# Mécanique

Statique

# *Mechanics*

*Static*

Réf :  
322 006

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
### Bending the pin

Shorten if necessary

**Knot to attach the string to the pin ('knotless knot'):** the reel of string is supplied.

Pass the string of the free snell through the eye of the hook.

Make a minimum of six rotations around the stem of the hook. The ideal is to cover the entire stem up to the curve of the hook.

Tighten

## 2 Additional accessories

- Magnetic table
- Box of marked weights with a hook

### 3 Operations

#### 3.1 Moment of a force in relation to an axis

##### 3.1.1 Concept of a moment

With the plate in any state of equilibrium, observe that by placing a very light object (for example, a single pin) in any of its perforations, the equilibrium is not modified (on account of inevitable friction) if the hole is close to the axis of rotation, whereas the plate rotates if the hole used is further from the centre.

Use varying lines of holes (horizontal, oblique, vertical) to demonstrate the action of a force, with rotation of a moving solid body around an axis  $\Delta$ , depends on the distance 'd' from its support to  $\Delta$ , and not the distance between the axis – point of application of the force.

Verify that irrespective of the intensity of the force, the latter does not modify equilibrium of the plate if the line of its effect meets the axis of rotation.

**Conclusion:** For a same applied force, the body rotates more or less depending on the position of the point of application of the force in relation to the axis of rotation. → The moment of the force is different.

##### 3.1.2 Expression of the moment


Apply ballast to the plate by applying a force  $\vec{F}_A$  to any pin (at A).

Using a second pin B, apply a different force  $\vec{F}_B$ .

Once the plate is immobile, mark the position of equilibrium by the values:

$\delta$  = distance between pin A and axis of rotation of the plate

d = distance between pin B and axis of rotation of the plate.

Demonstrate that if force  $\vec{F}_B$  retains the same intensity and the same direction, then movement of the pin around the axis of the force (at C or D on the above diagram) does not modify the equilibrium of the body.

Pin B can be placed in any other position and the intensity or direction of the force (or both) modified to restore the same equilibrium: observe that the product of  $F \times d$  are constant and characterise the action of a force, with equilibrium of the moving plate around the axis.

##### Expression of moments:

Moment at A:  $M_A = \delta \times F_A$  wh

where  $F_A$  is the force applied at point A and  $\delta$  the distance between the point of application of the force and the axis of rotation of the body.

Moment at B:  $M_B = d \times F_B$  where  $F_B$  is the force applied at Point B and D the distance between the point of application of the force and the axis of rotation of the body.

## 3.2 Theorem of moments

### 3.2.1 Addition of moments

- Retain the conditions for the previous experiment:
  - Same pin A, same force  $F_A$  (same weight).
  - Same equilibrium (same distance  $\delta$ ) as previously.
- Simultaneously apply two forces  $\vec{F}_1$  and  $\vec{F}_2$  using two pins  $B_1$  et  $B_2$ . Measure the distances  $d_1$  and  $d_2$ .
- Taking into account the signing convention, evaluate the moments:
$$M_1 = F_1 \times d_1$$
$$M_2 = F_2 \times d_2$$
- Compare the sum of  $M_1 + M_2$  to the value of the product  $M_B = F_B \times d$  in the previous experiment.
- Make several combinations of forces  $\vec{F}_1$  and  $\vec{F}_2$  corresponding to the same equilibrium of the plate, the moments  $M_1$  and  $M_2$  being of the same sign or opposite signs.

### 3.2.2 Condition for equilibrium of the body

Prepare an inventory of the forces applied to the body.

Verify that:

$$M_A = - M_B \rightarrow M_A + M_B = 0 \rightarrow M_A + M_1 + M_2 = 0.$$

Generalise the previous result and set out the theorem of moments:

If a solid body moving around an axis is in equilibrium, then the sum of the moments of forces which make it rotate in one direction is equal to the sum of the moments of forces which make it rotate in the other direction.

## 4 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**+33 (0)2 32 29 40 50**

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*Rubrique FAQ*

## Direct connection for technical support

**A team of experts at your disposal from Monday to Friday (opening hours)**

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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