Mécanique

Statique

Mechanics

Statics

Ref :
322 014

Français – p 1

English – p 9

Version : 9001

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
Apparatus set up with pulleys and overload.

Note the arbitrary orientation of the metal arms that carry the pulleys and the direct reading to within one mm on the slewable plate.

This apparatus is used to demonstrate and understand in a first approach to statics – both during courses and lab work – the role of the various parameters involved in the equilibrium of an object moving about its centre of gravity.

- Observation of moments
- Forces from masses or dynamometers

The forces F can be applied very easily at numerous points on the moving member. These forces can be :

- vertical (weight of the overloads);
- any direction, slewable by rotating an independent metal arm

The value of the distance “ d ” of the rotation axis to the support of the force is read directly and accurately on the graduated plate.

- Study of equilibrium
- Notion of couple
- Theorem of moments

Note : The experiments described herein can be performed either with masses or with dynamometers.


1 Principle – description

1. Moving member –scale –arm assembly, mounted on a nickel-plated rod (dia. 14 mm).

The moving member is a PVC plate with an arbitrary shape mounted on a highly sensitive ball bearing. It is perforated over its whole surface with holes 15 mm apart along two rectangular directions that can be used to locate the coordinates of a hole. The plate is ballasted so that its mass centre is on the rotation axis Δ (indifferent equilibrium).

A slewable rectangular metal plate with hydraulic friction carries two parallel graduations in mm (80, 0.80) on its smaller sides. The zero line passes at any time through the rotation axis Δ , *this scale allows the immediate measurement of d .*

Four nickel-plated metal arms separately mounted on a small square plate are independently slewable with hydraulic friction. The curved end of the arms has 2 holes in which a rod 6 mm in diameter can be inserted.

2. Pulleys, Ø 50 mm with highly sensitive ball bearings, mounted on a rod Ø 6 mm. In case of wobble, recouple by hand the PVC disk onto the bearing cage.

3. Dynamometers, 1 N and their clips mounted on a rod (Ø 6 mm). Place the mass centre of the dynamometer on the clip axis.

4. Very light studs (mass < 1 dg) to be inserted into the plate perforations and connected by means of wires to the plates (mass 20 g) that carry the overloads (20 g each) or to dynamometers.

Ensure that the wire that materialises the line of action of the force \vec{F} under study is perfectly parallel to :

- The plane of the plate, by driving the pulley rod more or less in ;
- The graduation lines of the plate oriented for that purpose during the measurement of the distance “d” of the line of action of \vec{F} relative to the axis Δ . Reduce parallax errors by driving the studs fully in.

1.1 Composition of the system P/N 322 014

- 1 moving member – scale – arm assembly
- 4 pulleys with 4 studs and 4 wires
- 4 dynamometers, 1 N with 4 clips, 4 studs and 4 wires.
- 4 overload carrying plates and 12 overloads of 20 g each.

Important : the ball bearings are used dry in order to reduce the friction forces. Do not lubricate them.

1.2 Additional accessories

- 1 universal stand (P/N 701 030 or 701 031) and 1 fastening fixture (P/N 703 245) to secure the apparatus spindle horizontally.
- 1 box of marked masses with hook (P/N 703 016) totalling 550 g.

1.3 Spare accessories

- Pulleys P/N 323 017
- Dynamometers P/N 322 049
- 6 sets of wires and studs P/N 323 021
- Arbitrary figure P/N 323 023
- Set of 4 dynamometer spindles P/N 323 022

2 Experiments

2.1 Moment of a force relative to an axis

2.1.1 The notion of moment

With the plate at indifferent equilibrium, insert a very light object (e.g. a small piece of match 15 mm in length) into one of the plate perforations and observe that the equilibrium is not modified (because of the inevitable friction) if the selected hole is close to the axis, and that the plate turns if the selected hole is farther off centre.

Use different rows of holes (horizontal, oblique, vertical) to demonstrate that the action of a force on the rotation of a solid moving about an axis Δ is dependent on the distance "d" of its support relative to Δ and not to the distance between the axis and the point of application of the force.

Verify that, irrespective of its intensity, the force does not alter the equilibrium of the plate if its line of action meets the axis $d = 0$.

2.1.2 Expression of moment

Ballast the plate (Apply a force \vec{F} to an arbitrary stud)

Example:


$\mu = 0,20 \text{ N}$;

$A(X = 5, Y = -3)$.

Using a second stud B apply a force \vec{F} . Once the plate is stationary, locate its equilibrium position by the value of δ . Measure the distance d.

If the direction of the force \vec{F} meets several holes B, C... demonstrate that, for the same intensity and direction of the force, moving the stud from B to C... will not alter the equilibrium (experiment easily performed with an oblique force).

Use other arbitrary positions for stud B and modify the intensity or direction of the force (or both) to restore the same equilibrium : demonstrate that the products $F \cdot d$ are constant and characterise the action of a force on the equilibrium of the plate moving about the axis.

2.2 Theorem of moments

2.2.1 Addition of moments

Ballasted plate: same stud A, same force \vec{f} ,
same equilibrium (same δ as above).

Simultaneously apply two forces \vec{F}_1 and \vec{F}_2 via two studs B₁ and B₂.
Measure the distances d₁ and d₂.

By following the sign convention, evaluate the moments:

$$M_1 = F_1 \cdot d_1 \quad \text{force } F_1 ;$$

$$M_2 = F_2 \cdot d_2 \quad \text{force } F_2 ;$$

Compare the sum $M_1 + M_2$ to the value of the product $M = F \cdot d$ obtained in the previous experiment.

Use several combinations of forces \vec{F}_1 and \vec{F}_2 corresponding to the same plate equilibrium, the moments M_1 and M_2 having identical or opposite signs.

2.2.2 Plate equilibrium condition

List the forces applied to the plate by taking force f into account. Calculate its moment $M' = f \cdot \delta$.

Verify that : $M' = -M$;
or : $M + M' = 0$;
or else $M' + M_1 + M_2 = 0$

Generalise the above result and express the theorem of moments.

In both cases the row of holes xx' is horizontal : same equilibrium position for the plate (same value of $\delta = 75$ mm)

2.3 Couple

For this study, ballast the plate and keep the same equilibrium: same stud, same mass μ , same value of δ .

2.3.1 Lever arm of the couple (figure 3)

Using a stud B_1 , apply a force \vec{F}_1 such that the support of \vec{F}_1 is parallel at equilibrium to one of the perforation rows, or YY' .

Select a hole B_2 in the row YY' passing through the rotation axis Δ , insert a stud therein and apply a force \vec{F}_2 parallel to \vec{F}_1 and having the same intensity but an opposite direction (ensure that the wires corresponding to \vec{F} and \vec{F}_2 are parallel to the graduations on the slewable plate).

The equilibrium is not modified since $d_2 = 0$. It is maintained when, with \vec{F}_1 and \vec{F}_2 remaining parallel and of the same intensity, the lever arm of the couple, D , remains constant, e.g.:

- With \vec{F}_1 and \vec{F}_2 kept parallel to YY' , move B_1 and B_2 by an equal number of holes in the direction XX' and in the same direction;
- If the common direction becomes parallel to XX' , maintain B_1 , the position of B_2 can be found easily.

2.3.2 Expression of the couple

Apply in succession the following couples with arbitrary directions :

F_1, F_2 lever arm D
 F'_1, F'_2 lever arm D'


Fig. 2


Verify that:

$$F_1.D = F'_1.D' = C$$


These products represent the respective moments of the applied couples since this moment is reduced to the moment of one of the forces when the support of the other force meets the axis.

2.3.3 Addition of couples

Apply two simultaneous couples and verify that the sum of their moments is C.

2.3.4 Generalisation of the theorem of moments

$$A \begin{cases} x = -4 \\ y = -3 \end{cases}$$

$$\begin{aligned} \mu &= 30 \text{ g} \\ \delta &= 63 \text{ mm} \\ m_1 &= m_2 = 42 \text{ g} \\ d_1 &= 45 \text{ mm} \quad d_2 = 0 \quad \{ \\ d_1 &= -15 \text{ mm} \quad d_2 = 30 \text{ mm} \quad \{ D = |d_1 - d_2| = 45 \text{ mm} \end{aligned}$$

Same plate equilibrium in both cases.

2.4 Measurement accuracy

The order of magnitude of the approximation is:

1 dg on the value of masses (stud, friction);

5 % on the value of the forces measured with dynamometers;

0,5 mm on the measurement of the lengths d, D

Hence an accuracy of approx. 2-3% on the value of the moment of a weight,
and approx. 6% in the case of a dynamometer.

3 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries