Mesure

Measurement

**Réf :
322 033**

Français – p 1

English – p 4

Version : 8006

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1 General

1.1 Purpose of the apparatus

The transparent labels, glued onto the four sides of a 10 N dynamometer (reference 322 050), allow measuring the weight of a body of mass m ($0 < m < 1\text{kg}$) on the Earth and to simulate (on the Earth) what its weight would be on the Moon, Mars or Jupiter.

If the mass of a body is constant from one planet to another whereas its weight varies, this means the force from a distance called weight and exercised by a planet on the body, depends on the planet itself, that is on its mass and size.

1.2 Description

The labels are presented :

- For the Earth: the name 'Earth' only, to be glued above the original graduations of the 10 N dynamometer reference 322 050.
- For the Moon, Mars or Jupiter: the name and the graduations corresponding to the weight of an object of a mass less than or equal to 1 kg.

	TERRE	LUNE	MARS	JUPITER
	0	0	0	0
	1	0,16	0,37	2,3
	2	0,32	0,74	4,6
	3	0,48	1,11	6,9
	4	0,64	1,48	9,2
	5	0,8	1,85	11,5
	6	0,96	2,22	13,8
	7	1,12	2,59	16,1
	8	1,28	2,96	18,4
	9	1,44	3,33	20,7
	10	1,6	3,7	23

Original dynamometer graduations

1.3 Technical characteristics

Transparent self-adhesive pre-cut labels (104 x 13 mm) adapted for the 10 N dynamometer ref. 322 050.

2 Instructions for use

2.1 Precautions for use

The labels are self-adhesive. Since their glue is permanent, check the 'zero' graduations of each planet are at the same level as the original 'zero' on the dynamometer.

2.2 Use

Attach a marked mass to the end of the dynamometer. Measure the weight using the original graduations and by simulation, read off the weight the marked mass would have if it was on the Moon, Mars or Jupiter.

Remark : This is indeed a simulation for if it was possible to go to Mars with a dynamometer, for example, for a mass of 1 kg, the weight 3.7 N would be read on the original graduations of the dynamometer and not on the transparent labels corresponding to the planet 'Mars'.

Taking into account the proportionality relation between weight and mass: $P = m.g$, it is possible to calculate the intensity of the gravity g of each planet.

3 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries