Mécanique

Dynamique

Mechanics

Dynamics

Réf :
332 052

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Ballistic Launcher

Français – p 1

English – p 9

Version : 5108

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1. Description

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0 825 563 563*

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Fig. 12- 3D image of the ballistic launcher

1. Description

1.1 Photograph of the ballistic launcher

Fig. 2- Ballistic launcher fixed to a bench with its bench mounting

1.2 Contents of the packaging

- The launcher with inclination gauge
- Bench mounting stem
- Clamp for bench mounting
- Accessory for horizontal propulsion of a carriage (or something else)
- Propulsion accessory for a ball-bearing (supplied: fitted in the launcher)
- 3 ball-bearings


Fig. 3 Components of part number 332052.

To use the carriage propulsion accessory:

The pictures below illustrate the operations to be performed.

Fig. 4 Remove the two rubbers caps from the pin linking the two springs. Slide the pin out of the launcher.


Fig. 5 Remove the ball-bearing propulsion accessory. Fit the carriage propulsion accessory by first inserting the metal part.


Fig. 6 Through the side openings in the launcher, insert the spring-holding pin passing through the orifice in the accessory. Refit the rubber caps.


Fig. 7 Fit the springs so as the hold this rod.

1.3 Additional equipment

1.3.1 Essential equipment

1.3.1.1 Power supply

For any manipulation of the ballistic launcher, the experimenter will require a minimum power supply of 12 V – 0.2 A.

EVO2® F6F12 power supply: 281483; or EVO2® F3F12 power supply part number 281404.

Operating procedure:

Attach and orientate the ballistic launcher.

Place a ball-bearing in the launcher.

Connect the ballistic launcher to the power supply.

1.3.1.2 Set the power supply to 12 V.

Open the circuit to trigger firing.

Warning: do not place fragile objects in the line of fire.


1.3.1.3 Video recording device

To record any movement of a projectile sent with the ballistic launcher, we recommend using a fast camera.

1.3.2 Necessary equipment

1.3.2.1 Mechanical bench

The ballistic launcher is supplied in the 'mechanical bench' pack, part number 332051. Experiments connected with this bench are shown in the instructions for the bench.

*Fig. 8- Ballistic launcher and free fall mounted on the mechanical bench
The monkey and the hunter experiment*

1.3.2.2 Free fall

With the free fall – part number 332023 and a Modumontage® [retort stand] device, carry out an experiment to demonstrate conservation of the quantity of movement.

Commonly called the 'monkey and hunter' experiment, it consists of:

- Fit the free fall
- Attach the ballistic launcher
- Mount in series: the free fall, launcher and power supply
- Aim the ball-bearing from the free fall with the ball-bearing from the ballistic launcher

The difficulties with this apparatus are successfully aligning the two devices (see photos fig. 9). To facilitate this experiment, we recommend setting up the apparatus shown in the picture below.

- Attach a retort stand rod to the 2 feet of the free fall
- Attach the launcher to the rod


Fig. 9- Assembly of the Free Fall + Ballistic Launcher for the 'monkey and hunter' experiment.


To assemble this apparatus, the experimenter will require the following additional equipment:

- 1 Modumontage[®] rod, L700 mm – part number 703459
- 1 Modumontage[®] clamp – part number 703529

2. Examples of experiments

2.1 Ballistic firing for measuring gravity

The simplest experiment with this apparatus is to perform a ballistic firing. To do so:

- select the inclination of the ballistic launcher using the scale on the launcher
- connect the launcher to a 1 V power supply (e.g. part numbers 281483 or 281404)
- place a camera (fast camera type: 572 001, with software) so as to view the trajectory of the ball-bearing

Example of trajectory obtained:


Fig. 10- Trajectory obtained with Jeulin ballistic launcher and fast camera + Cinéris software (reference 572001)

The graph shown in figure 5 is a parabola. Its equation, given the launcher's angle of inclination, is used to find the value of the gravitational field.

2.2 Conservation of the quantity of movement

The results of the 'monkey and hunter' experiment are given in figure 6.


Fig. 11- Result by chronophotography of the 'monkey and hunter' experiment. For 3 firings at different velocities. The red ball-bearing comes from the ballistic launcher (on the right of the picture). The silver ball-bearing comes from the free fall.

This experiment shows that after the impact, the trajectories are reversed. In this way the conservation of the quantity of movement is demonstrated for 3 different initial firing velocities.

2.3 Measurement of flight time

Using the impact detector – part number. 332024 and a counter, it is easy to determine the flight time of a ball-bearing.

Description of the apparatus:

- Attach the ballistic launcher to a bench using the clamp provided.
- Place the ball-bearing in the launcher and connect an appropriate power supply to the launcher.
- Place the impact detector at the position of the ball-bearing's point of impact.
- Connect the launcher and impact detector to the counter.
 - o On the launcher, this is an opening
 - o On the detector, this is an closure

2.4 Precautions for use

Since the firing velocity exiting the launcher is quite high (given experiments on the mechanical bench), it is essential to take precautions before turning on the air.

For example, avoid standing in front of the cannon when firing. To do so, always place the ball-bearing in the launcher after adjusting all the settings. In this way, you will avoid any accident.

In the same way, it is highly recommended not to leave anything on the bench when firing (telephone, laptop computer, tablet, lamp, etc.).

3. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux