Mécanique

Dynamique

Mechanics

Dynamics

Réf :
332 072

Banc de mécanique complet 1m20
Complete 1.20 m mechanics bench

Français – p 1

English – p 10

Version : 6106

1 Description

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0 825 563 563*

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1 Description

This complete and moduable assembly can be used to perform experiments on a wide variety of translation and parabolic movements. The rigid aluminium bench is graduated and equipped so that it can be placed either horizontally with a level adjustment, or inclined with an additional stand. Three sights are provided to avoid parallax errors during a video recording.

Many experiments are carried out comfortably and easily:

- Carriage movement by different propulsion modes (spring, constant force, manual action).
- Impact of two carriages (conservation of movement quantity).
- Determination of vectors: speed and acceleration during a uniform or accelerated movement of a carriage.
- Determination of mechanical, kinetic and potential energies.
- Ballistic by means of the supplied launcher (parabolic movement)
- Relativity of the movement, study of the movement in a coordinate system.

2 Composition

The mechanics bench is composed of:

- 1 graduated bench length 1. 20 m. Adjustable stands and stops
- 2 secured carriages to prevent them from dropping
- One of the two carriages is fitted with an exploder
- 1 ballistic launcher ref. 332 052
- 1 table attachment device for launcher
- 1 set of accessories for attachment of forks and free fall rod
- 1 pulley on support
- 1 string for pulley
- 1 mass (about 500g) to be fixed on the carriage
- 2 rods cross-section 15x15 mm one of which is silk-screen printed + fasteners
- 3 Modumontage[®] vice clamps
- 3 sights


Sight

Nut attachment to sight

Vice clamp to fix launcher to table (or working surface)

Ballistic launcher

Mobile with rod

Ballistic launcher on bench

3 Usage

3.1 Installing the bench

The bench shall be placed on a table with a uniform background colour to give good contrast during a video recording.

It shall be checked that the bench is horizontal by placing a spirit level on it, or by putting the carriage on it and checking that it doesn't move.

The adjustment is made by screwing / unscrewing support legs.

A stand (not supplied) should be used with the support rod sleeve accessory (supplied) for manipulations on an inclined plane. See section 3.6.3.

Adjustable legs

A coloured sticker shall be fixed on the side of the carriage and is used as a marker for video positioning of the carriage. The three sights (magenta coloured crosses) must then be positioned in the vertical plane containing the coloured sticker. This precaution is essential to avoid parallax errors due to the webcam and to obtain precise measurements of the coordinates.

The sights are positioned along the two directions in the plane (vertical et horizontal, defining distance scales and the coordinate system of the plane). They will be used for positioning when processing the videos. The measurement plane may be different from the bench, depending on the experiment being carried out.

Carriage stops are mounted by clamping two clips on each side of the bench (see photos below). End pieces (three maximum) can be fixed on a stop depending on the experiment to be done.

Adjustable leg and sight

3.2 Ballistic launcher

The launcher is composed of the following elements:

- Exploder
- Inclinometer
- Attachment system (attachment rod to working surface and clamp)
- Power supply cables to start
- Ball propulsion accessory
- Carriage horizontal propulsion accessory
- 3 steel balls: $\varnothing 16\text{mm}$

The launcher can be used either fixed to a bench or in carriage propulsion mode. It can also be fixed to a table to launch a projectile (ball), or mounted on a carriage that may be immobile or moving.

Arm the launcher

A power supply must be used for all manipulations with the launcher: 12 V – 0.2 A minimum.

After fixing and then orienting the launcher, place a ball in the launcher and then switch the power supply on, and open the electric circuit to trigger firing. Either the ball propulsion or carriage propulsion accessory will be used depending on the experiment.

Note: the thrust is adjustable.

See section 3.3.1 for further information.

3.3 Mobiles

Carriages are made using the same materials as the bench. They are fitted with four wheels mounted on high quality ball bearings to minimise friction. The carriages fit the bench rail perfectly with no lateral friction.

These small low-mass ball bearings have negligible inertia compared with that of the moving carriage.

A coloured sticker to be glued on the side of the carriage is used to identify carriage positions and to make abscissa measurements on the images of videos made during the manipulations.

Carriage characteristics:

- Mass of carriage alone: about 500g
- Mass of carriage with exploder: about 530g
- Additional mass: about 500g
- Carriage length: 160 mm

For better precision, measure the mass of your mobiles with your equipment: additional mass, launcher or other.

Adjust the spring compression of the carriage fitted with the exploder. The spring compression can be adjusted using a screw on the rod.

There are two possible positions. Triggering is controlled using the push button on top of the carriage.

Compression adjustment screw

Two adjustment notches and trigger button

3.3.1 Propulsion

The following images illustrate operations to be done to install the propulsion accessory on a carriage:


Remove the two rubber end pieces from the pin connecting the two springs. Slide the pin outside the launcher

Take out the ball propulsion accessory, then put the carriage propulsion accessory into place, inserting the metal part first

Insert the spring retaining axis through the lateral openings of the launcher, passing through the accessory orifice. Put the rubber end pieces back into place

Put the springs into place to retain this rod

3.3.2 Set up the ballistic launcher


Launcher with recovery device

Launcher attachment screw

The ballistic launcher is fixed to the central groove of the carriage by a screw.

Different uses of the ballistic launcher

3.3.3 Propulsion integrated into the carriage

Adjustment of the propulsion control spring

There is a screw for adjusting the position of the spring depending on use.

3.4 Pulley


Pulley on support

Nut for attachment under the bench

The pulley must be installed on a support fixed to the central groove underneath the bench, and an optical barrier (not supplied) may be fixed to the same support.

3.5 Carriage magnet support piece


Carriage magnet support piece

Magnets can be removed from their support so that either attraction or repulsion with the second carriage is possible.

3.6 Accessories

3.6.1 End piece with magnet support

A magnet support piece can be fitted in the groove on the bench and act as a limit stop. Magnets can be removed from their support such that there is either an attraction or a repulsion force, to demonstrate movement of a carriage under the action of a magnetic force.

End piece with magnet support

3.6.2 Support for light barrier

Light barrier supports (not supplied) are fixed in the bench lateral groove.

Light barrier support

3.6.3 Bench inclination


Bench inclination accessory

There is a sleeve that fixes the bench onto a support stand to incline it.

3.6.4 Additional mass


An additional mass of about 500g can be installed on the carriage in its groove using the screw supplied.

Attachment of the mass on the carriage

3.6.5 Installation of the rod and optical forks

Put the optical forks onto the rod and then put the rod in the installed support. Use the Jeulin forks with the chronocounter or use the Jeulin forks in CAEX with the Foxy®, Tooxy® and AirNeXT® consoles and the chronoméca module made by Atelier Scientifique.

Installing the fork on the bench

4 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux