Mécanique

Mechanics

Réf :
332 127

Français – p 2


English – p 11

Version : 4211

Chute libre

Free fall

1 Description

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
1 Description

This device can be used to study linear motion under uniform acceleration. By considering the case of free fall under gravity, students will determine the acceleration, which, in this simple case is known and equal to the value of the earth's gravitational field: g (in ms^{-2}).

Compatible with our datalogging range, or any digital stop watch, experiments for determining " g " are described in section 2.3.

Compatible with our mechanics range, this product will allow you to demonstrate the conservation of mass as part of a collision experiment using our standard equipment, ballistic launcher and / or mechanics experiment bench. This experiment is described in section 2.3.5.

2 Apparatus and assembly

2.1 Apparatus

The product comprises (see image above):

- **1 x FREE FALL TRIGGER**

- 1 free fall trigger with 2 ball-bearing retaining tongues of which 1 is a spare. The ball-bearing must be positioned so that it is in contact with the 3 screw heads.
- 1 battery for ensure synchronization with ExAO sensors

- **1 x COLLISION DETECTOR**

- 1 collision detector with banana jacks

- **1 x ATTACHMENT DEVICE**

- 1 rod of cross-section 15x15 mm and length: 40cm
- 2 modular assembly feet
- 4 knurled screws to secure the modular assembly feet
- 2 modular assembly vice grips
- 1 rod of cross-section 15x15 mm and length 1m

Note that this 1 m rod is printed with two graduated rulers.

One of the printed rulers has its zero at the top of the rod. It enables work in Ex.A.O with the dedicated **chronoméca** software.


The other printed scale has its zero at the bottom. This zero corresponds exactly with the plane of contact of the ball-bearing with the collision detector.

- **2 x BALL-BEARINGS \varnothing 16mm**

2.2 Step by step assembly and reading of the positions

Referring to the image below:


- Secure 2 knurled screws beneath each modular assembly foot.
- Secure the 40 cm rod to the 1 m rod using the modular assembly vice grip



Ensure that the 1m rod is positioned at the floor without resting on it.

- Position the collision detector in front of the rod. Thus this zero corresponds exactly with the plane of contact of the ball-bearing with the collision detector.
- Fix the free fall trigger on the rod: 2 possible cases:

Experiment without datalogging system:


Experiment via a datalogging system with optical forks:

During experiments using optical forks, it is essential to have optical forks (453026) and modular assembly vice grips (703529) available:

Reading the fork
position

Reading the position: The green arrows on the images above show the best way of reading the position of the ball-bearing at each step of its trajectory.


The error in the distance travelled by the ball-bearing can be affected by a significant angle between the vertical and the rod. To avoid this error, the use of a plumb line, spirit level or inclinometer is recommended.

2.3 Use of the free fall trigger

The free fall trigger, is a mechanical device that requires a battery (1.5 V, included), more specially when using voltmeter sensor.

2.3.1 Triggering of the fall of a ball-bearing

The triggering of the fall of the ball-bearing is implemented by pressing the push-button as shown below:

The ball-bearing must be positioned in its slot to enable the correct operation of the trigger.

2.3.2 Measurement with a chronometer

Equipment required:

- Free fall system ref. 332127
- Chronometer ref. 351028
- 12V power supply unit ref. 281512

To perform measurements using the chronometer:

- Connect the free fall trigger on the left part of the chronometer.
- Trigger the start of the fall over an "opening". (See chronometer).
- Connect the collision detector to the right part of the chronometer. (See chronometer).
- Trigger the end of the fall over an opening.

The time displayed by the chronometer screen corresponds to the time necessary to traverse the distance separating the trigger and the receiver.

Based on the equation: $y = \frac{1}{2} \times g \times t^2$, determine the value of the acceleration due to gravity g .

By reproducing this experiment for different heights, plot the change in height as a function of the square of the time. From this determine the coefficient $\frac{1}{2} \times g$. Determine the value g .

2.3.3 Measurement using a data acquisition console

Equipment required:

- Free fall system ref. 332127
- 1 data acquisition console with synchronization (Foxy, AirNeXT, etc.)
- Optical fork ref. 453026
- modular assembly vice grips ref. 703529 (1 per fork)

To perform data acquisitions with the data acquisition console:


Without external synchronisation:

- Connect the free fall trigger to a voltmeter sensor on the input 1 of your console
- Connect the collision detector to input 2 of your console
- Synchronise your acquisitions at input 1, on a falling edge for the value less than voltage of the battery (approx.. 1.5 V).

With external synchronisation:

- Connect the free fall trigger to the external synchronisation element of your console
- Connect the collision detector to input 2 of your console
- Synchronise your acquisitions at the external synchronization device

Example of the signal acquired:

Read-off the delay elapsed between the start of the data acquisition and the first marked peak.

Based on the equation: $y = \frac{1}{2} \times g \times t^2$, determine the value of the acceleration due to gravity g .

By reproducing this experiment for different heights, plot the change in height as function of the square of the time. From this determine the coefficient $\frac{1}{2} \times g$. Determine the value g .


2.3.4 Data acquisition with a data acquisition console and optical forks

Equipment required:


- 1 Free fall system: 332 127
- 1 data acquisition console Foxy®, AirNeXT®, Tooxy® etc.
- 1 to 3 time-synchronized forks: 453026
- 1 to 2 time-synchronized adapters: 482048
- 1 DIN banana connector lead for the free fall system: 283628

The experiment shown here was performed using one time-synchronized fork and 1 time-synchronised adapter.

Photographs of the assemblies of the free fall trigger and 1 optical fork:

The results obtained are presented here by displacement of the optical fork along the whole length of the trajectory.
Note the collision detector can be used to receive the ball-bearing, but it cannot be used for determining the limit of travel.

Example of the experimental result using a fork and chronoméca software

In the left insert, slide the optical fork(s) on the vertical axis of the free fall device. Enter the height value for each fork.

The result obtained in figure 3 was obtained by moving the fork along the vertical axis. Successive data acquisition points make it possible to view the parabola representing the curve: $y = \frac{1}{2} \times g \times t^2$. The coefficient of this parabola yields the value of the gravitational field.

N.B.: The experiment described in this section is also compatible with AirNeXT® in connected mode.

2.3.5 Experiment: elastic collision

2.3.5.1 Reminder

Elastic collision:

During an elastic collision between 2 solid, non-deformable ball-bearings, momentum is conserved: this is referred to as an elastic collision.

Theory of the conservation of momentum:

Based on Newton's second law $\sum \vec{F} = m \cdot \vec{a} = m \cdot \frac{d\vec{p}}{dt}$

In the case of a system where there are no net external forces or the forces cancel each other out:

$$\sum \vec{F} = \vec{0} = m \cdot \vec{a} = m \cdot \frac{d\vec{p}}{dt}$$

From which is deduced: $\frac{d\vec{p}}{dt} = \vec{0}$.

The momentum (represented by p) does not vary over time, i.e. it is conserved.

$$p = p_1 + p_2 = \text{constant.}$$

2.3.5.1 Experiment

Equipment required:

- Free fall system: 332 127
- Ballistic launcher: 332052

Or Mechanics experiment bench: 332072

- 1 12V DC power supply: 281483
- 1 high-speed camera: 572000

Supplementary equipment recommended to facilitate the alignment:

- 1 rod of cross-section 15x15 mm and length 70cm
- 1 modular assembly vice grip

Assembly:

Secure the feet of the free fall system, position the 70 cm rod horizontally. Using a modular assembly vice grip, secure the ballistic launcher as shown in the photograph below.

Photograph of the complete assembly



Photograph of the ballistic launcher assembly

Connections:

Connect the launcher, the free fall system and the power supply in series. Thus when the fall of the ball-bearing is triggered, the launcher is likewise triggered. The two events are synchronised.

It is thus possible to demonstrate the conservation of momentum by varying the firing force of the launcher.

3 firing speeds were processed. For each speed, a video was recorded using the high-speed camera (reference: 572000). The chronological photographs shown below result from these captured videos. The software used is Cinéris.

Chronophotography of the elastic collision experiment for shots with 3 different speeds.

These experimental results show that:

- Whatever the speed of firing the two ball-bearings meet
- The trajectories of the two ball-bearings are inverted.
- The speed of the red ball after the collision equals the speed of the grey ball before the collision. Conversely, the speed of the grey ball after the collision equals the speed of the red ball before the collision.

This experiment is also called: the experiment of the monkey and the hunter.

N.B.: To be able to clearly visualise the inversion of the trajectories, it is recommended that one of the two ball-bearings is coloured.

To extend the experiment, it is possible to use balls of different masses. Ensure that the ball-bearing used in the free fall trigger is conductive.

3 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops. For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 23