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Periodic phenomena

**Ref :
351 043**

Français – p 1

English – p 3

Version : 8009

Stroboscope tachymètre digital

Digital tachometer stroboscope

1 Principe – description

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1 Principe - description

This stroboscope is a source of light with brightness particularly suitable for the observation and measurement of periodic motion.

The device is equipped with four control buttons:

- One ON-OFF switch,
- One range selector,
- One flash frequency control button,
- One fine adjustment button for the frequency of flashes,
- One display (4 digits) of the number of flashes/min,
- One 2 P + T power outlet.

2 Technical characteristics

- 220 V, 50 Hz \pm 10 % Power supply,
- Xenon lamp 40 joules,
- Frequency range of flashes,
 - 100 to 1000 flashes/min,
 - 1000 to 10 000 flashes/min.
- Display accuracy: \pm 1 flash from 100 to 5000 flashes/min,
- Display accuracy: \pm 0.05 % more than 5000 flashes/min,
- Dimensions: 210 x 120 x 120 mm,
- Mass: 1 kg approximately.

3 Commissioning

- Connect the power cable to the stroboscope and to the 220 V mains,
- Turn "ON" the switch,
- Turn the device towards movement to be immobilised,

- Slowly turn the frequency button from the maximum position until a fixed image appears.

As soon as the first fixed and single image is obtained, read the frequency of the motion observed on the display.

From the fixed image, a slight modification of the position of the fine adjustment button in either direction can lead to a dynamic scanning by creating a slow apparent movement.

Note:

Operating time:

- below 3000 flashes/min: 30 min
- above 3500 flashes/min: 5 min

For an operation time greater than those indicated above, we advise you to keep a cooling period (about 10 minutes) so as to avoid damaging the device through over heating.

4 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
0825 563 563

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries