Mesure

Initio® 2

Measurement

Initio® 2

Réf :
351 071

Français – p 1

English – p 7

Version : 7105

Sonomètre Initio® 2
Initio® 2 Sound Level Meter

1. Description

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0 825 563 563*

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1. Description

The Initio® 2 sound level meter is a device that measures the sound level picked up on a microphone. It also measures the sound pressure and is equipped with a sound analyser. It can also be used to record data, store, and then restore it later on a computer.

A WI-FI module (ref 451 470) is also available for viewing and recording measurements on a tablet or a PC equipped with WI-FI.

Equipped with an analogue output, the Initio®2 sound level meter can be used with an ExAO interface equipped with direct inputs via banana jacks or via a Groomy® interface and jack.

2. Utilisation

Upon first use, it is advisable to completely charge the Initio® 2 battery. Use the μ USB (7) supplied in separate packaging to charge the battery. In general, we recommend using the instrument in complete charge-discharge cycles (wait until it is completely discharged before fully recharging).

2.1 Start-up

- Press the start button (1) situated on the left side of the instrument.

The measurement displays directly on the instrument's screen.

- To switch off the Initio® 2, briefly press the start button (1) and confirm by pressing the knurled knob (3).

2.2 Browsing

Initio® 2 has three navigation buttons:

Selection-confirmation (3):

The central knurled knob (3) is to scroll through the menus of the Initio® 2. Briefly pressing the "OK" button in the center of the knurled knob confirms the selection.

Back (2):

The "Back" button (2) is used to switch from one type of display to another from the main screen in the following order:

- Measurement
- Table
- Curve

Within a menu it is used to go back a step.

Menu (4):

This button displays Initio2®'s main menu.

2.3 Tacking measurements

Sound level meter mode

This mode has three measurable parameters:

- Noise levels in dB
- Weighted noise levels in dB (A)
- The effective sound pressure in Pa

To change the parameter, press the knurled knob from the main menu.

Analyser mode

This mode provides information on the sound frequency.

To access:

- Press the **Menu** button
- Using the knurled knob, select **Display** and press on the knurled knob.
- Select **Analyser** and confirm.

The information provided is:

- Above, the shape of the sound wave with the fundamental frequency calculated by detecting the "zero" crossings.
- Below is the sound spectrum (Fourier Transform) with the fundamental frequency obtained by taking the maximum intensity on the Fourier transform.
- Bottom left, the musical note corresponding to the calculated frequency.

To return to sound level meter mode:

- Press the **Menu** button
- Select **Display** and confirm with the knurled knob
- Select **Display** and confirm

2.4 Data acquisition

Initio® 2 has an internal memory that enables it to save measurements. Two modes are available:

- Manual acquisition:

Acquisition is on a point by point basis. Each measurement is manually recorded by the user (see below)

- Time acquisition:

Data acquisition takes place as a function of time, according to a duration entered by the user.

2.4.1 Data acquisition

- Press the Menu button (4)
- Using the knurled knob (3), select **Acquisition** and confirm by a short press on the knurled knob (3).
- Select **Manual**
- To record a measurement, select the **Save** icon (white circle) and confirm by briefly pressing the knurled knob (3).
- To stop the recording select the **Stop** icon (black square) and confirm by pressing the knurled knob (3). Confirm the selection in the dialogue box.

2.4.2 Time acquisition

- Press the Menu button (4)
- Using the knurled knob (3), select **Acquisition** and confirm by briefly pressing on the knurled knob (3)
- Select **Time**
- Using the knurled knob (3), select the number defining the duration and press on the knurled knob (3).
- Modify this duration by turning the knurled knob (3) then confirm by pressing on the centre of the knurled knob.
- If the data acquisition is to loop (return of the acquisition to zero at the end of a defined period) select **Loop**
- Confirm the final choice by selecting **OK**

Data acquisitions based on time last between 60 seconds and 50 hours. The interval between two measurements adjusts automatically.

2.4.3 Recording the data

Once the acquisition is complete (time or manual) it must be saved to be able to retrieve the points on a computer.

- Press the Menu button (4)
- Select **File** and confirm
- Select **Save**

Only a single data acquisition can be saved. By selecting **Open**, you open the last saved data acquisition.

2.4.4 Transferring the data to the computer

Once the acquisition is complete, it is possible to transfer the data to a computer to process it:

- Install the Initio® 2 software available at www.jeulin.fr/Initio2
- Connect the Initio® 2 to the computer via the USB-μUSB cable
- Launch the Initio® 2 software
- Load the saved acquisition by clicking on the icon:

The data then appears on the screen. Right-clicking the chart provides access to various graphic tools.

3. Data visualisation

It is possible to view the data at any moment in tabular or chart format. To do so, simply press the Back button (2). The display then toggles between table and chart measurements. Regardless of the current display, data acquisition can be monitored.

After 5 seconds without activity, the table or chart is displayed in full screen. To again obtain the associated functions, turn the knurled knob.

In the table, the arrows situated at the bottom right are used to move between the rows of the table.

4. Analogue outputs

Initio® 2 is equipped with banana jack and 3.5 mm jack analogue outputs.

4.1 Banana jacks


- Connect the Initio2 to an ExAo interface equipped with direct inputs using banana leads.
- Configure the direct input in the following manner:
 - o 0dB corresponds to -5V
 - o 120dB corresponds to +5V

Example with the Foxy® console

- Start the generalist

- Place the icon of a direct input on the desired axis
- Select the calibre $\pm 5V$

- Fill out the customisation tab as below:

- Click on the disk to save the configuration

Configuration is complete, the value displayed is now the sound level with symbol L and the units are dB.

4.2 Jack

The 3.5mm jack is to connect the Initio2 to a Groomy® interface (ref. 181 000) for example.
The voltage supplied varies between 0 and 5V. So that:

- 0 dB corresponds to 0V
- 120dB corresponds to 5V

5. WI-FI communication

Thanks to the WI-FI module (ref. 451 470) it is possible to view the Initio® 2 data on a tablet for example.

5.1 Installation of the WI-FI module

The WI-FI module is implemented with the Initio® 2 switched off.

- Raise the handle of the Initio2
- Remove the back by sliding it down
- Position the WI-FI module at the back
- Push the module firmly into place until it clicks. The handle must be able to retract without forcing.

5.2 Viewing on a tablet

Once the WIFI module is installed, the tablet must be configured:

- Press the Menu button
- Select Settings
- Select Network

The module information is displayed on the screen.

- On the tablet, choose the name of the WI-FI network displayed on the Initio® 2.
- Open a web browser
- Enter the address displayed by the Initio2 in the address bar of the browser.

6. Settings

Various information concerning the Initio can be obtained by pressing on Menu and selecting Settings → Information

This page displays the serial number of the device and the software version, the battery charge level and the operating language. The latter is changed by turning the knurled knob.

7. Update

On our website at the address <http://en.jeulin.fr/en-resources/english-resources/software/free-software.html>, the last update for your Initio® 2 is available. Its version is indicated in the designation.

In order to know the version of your Initio® 2, press the “Menu” button and then “Setup” “Informations”.

8. Technical specifications

Sound level meter mode

Sound level

- Range: 0 to 120 dB

Effective sound pressure

- Range: 0 to 20 Pa

Analyser mode

- Frequency range: 50 Hz to 16 kHz

Analog Output

-5 V to +5 V via 4 mm safety banana jacks

0 to +5 V via 3.5 mm jack

Communication

USB, WI-FI optional

Power

Li-ion battery

Dimensions

95 x 123 x 44 mm

9. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT

468 rue Jacques Monod

CS 21900

27019 EVREUX CEDEX FRANCE

+33 (0)2 32 29 40 50