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***Measurement and
instrumentation***

**Ref :
353 026**

Français – p 1

English – p 5

Version : 6010

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NOTES

1 Principle – description

When connected to an electronic timer, the time counter sensors allow time measurements :

- A first sensor starts the time measurement when an opaque object interrupts the beam ;
- The second sensor stops the time measurement when its beam is interrupted

Note : *In free fall experiments, only one sensor can be used. In this case, the timing will start e.g. upon opening of the circuit closed by the ball.*


These time counter sensors feature the advantages derived from the research on the computer-aided scientific experimentation system (ESAO) and they are high-performance devices :

- Very accurate thanks to the narrow beam and the “clean” electronic signal, they can be used in conjunction with 1/1000th-of second timers ;
- Sturdy and reliable, they are practically unaffected by ambient light as emission and reception are protected by a slit and infrared beams are used.

2 Technical data

- Banana plugs Ø. 4 mm
 - o Red : power supply hot wire
 - o Black : ground
 - o Yellow : signal - 0 Volt with uninterrupted beam
- Narrow IR beam
- Insensitive to ambient light
- Perfectly “clean” signal
- Suitable for the various electronic timers used in your laboratory

- IR beam interruption indicator
- Supplied with a 100-mm long dia. 6 mm threaded-end rod

- The output signal fed into the timer is 0 V when the circuit is closed, and 12 V when the circuit is open (red beam). The output signal is available on an open-collector transistor. A load of $470\ \Omega$ to $10\ \text{k}\Omega$ is required in the timer between the signal input and the power supply + pole. This allows you to adapt the time counter sensors to various timer models. JEULIN's time counter P/N 351030 contains a resistance within the above range. Therefore the sensors can be directly connected to this timer.
- Supply voltage : 4 to 6 Volts d.c. (3.5 V rms with the time counter P/N 351030).
- Current consumption : approx 100 mA.
- Connecting cable, 2 m, terminated by 3 banana plugs : red, black, and yellow.
- Dimensions : 110 x 65 x 15 mm (without the cable).
- Mass : approx 160 g

3 Experiments

These time counter sensors used in conjunction with an electronic timer allow accurate measurements of the travel time of an object (time, average velocity). Counting is also possible.

Experiment example :

1. Prepare the experimental setup : free fall, air cushion table, pendulum,...
2. Screw the threaded-end rods onto the sensors by hand (do not use a wrench). Generally two sensors are used with one electronic timer. Place the two time counter sensor in the path of the moving object. The accuracy of the results obtained is dependent on the sensor alignment relative to the moving object path. In free fall experiments, use a plumb line to align the ball starting axis and the sensor beam marks. When using an air cushion table, add opaque indexes (e.g. 1.5 mm in diameter) on the carriages(s) as required. Then position the sensors parallel to the table.
3. Connect the two time counter sensors onto the electronic timer.

Note 1 : If the power supply (+ and – poles) is not provided via the timer, connect the supply ground (4 to 6 V d.c.) to the timer supply ground.

Note 2 : *When using an electronic timer not supplied by JEULIN, measure the voltage between ground and input. If this voltage is approx. 0, add a $470\ \Omega$ to $10\ k\Omega$ resistor between the power supply + pole and the signal input inside the timer box (add approx. $470\ \Omega$ for a TTL input).*

4. The indicator on the sensor side comes on when the beam is interrupted. This is used to verify that the beam is actually interrupted by the moving object.

The first counter sensor can also be placed exactly at the starting point of the moving object : with the moving object attached to its starting point (ball, carriage index, etc.), bring the sensor near until the indicator comes on. Then move the sensor away and position it exactly at the point where the indicator come off. Then the sensor is located at the starting zero point ($t = 0$, $v = 0$, $L = 0$).

4 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*