**Ex.A.O.**

Console

***Data logging***

Console

Réf :  
451 461


Français – p 1

English – p 9

Version : 7112

**Console ESAO Visio<sup>®</sup> Plus**  
***ESAO<sup>®</sup> Visio Plus console***

## 1. Introduction

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
**0 825 563 563\***

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## 1. Introduction

ESAO® Visio Plus is a Computer Assisted Experimentation console using USB technology for increased speed and simplicity.

A specific feature is the option to connect a camera for synchronised data and image capture.

A Webcam or a fast Jeulin camera can be connected to the rear.

Complemented by a wide range of ESAO® sensors and adaptors, ESAO® Visio Plus means you can measure all types of magnitudes and perform a variety of experiments, notably thanks to the scientific workshop software collection, Science studio.

## 2. USB technology

This technology provides power comparable with a PCI; card but also and above all, unprecedented flexibility of use.

Compatible with most equipment, the ESAO® Visio Plus console means you can take advantage of all the benefits of USB technology (speed, ease of connection, etc.) so ensuring your equipment will be compatible with new developments.

## 3. Principle of capture


A sensor converts a physical magnitude into an electrical signal. This signal is then converted by the appropriate adaptor to make it compatible with the analogue input of the console.

The console performs analogue digital conversions, memorises the measuring points and sends these to a computer via the USB connection.

The computer then processes the measuring points and, for example, displays the results as a graph.

The console can also measure time and record events.

## 4. Description of the console

The console incorporates the various input/outputs accessible to the user. Some are accessible on sockets, others on connectors which can also accept adaptors.

### 4.1 Locations for adaptors

On the front panel of the unit, four locations can independently accept four adaptors from the entire ESAO® 4 range.

**On channels 1 and 2, two connectors can be used to connect two adaptors from the ESAO® 3 and ESAO® 3.1. range.**

### 4.2 Direct analogue inputs

Four  $\pm 5$  V differential analogue inputs are accessible on safety sockets (diameter 4 mm): a red one for the + pole and a black one for the - pole. These inputs are protected against permanent over-voltage to 250 V.

### 4.3 Analogue outputs

S1 designates the 0 to 10 V analogue outputs. It is accessible on the top of the console on a safety socket (diameter 4 mm) which is yellow. This output can supply DC voltage or act as a generator for periodic functions.

S2 is the  $\pm 5$  V analogue output. It is accessible on the top of the console on a safety socket (diameter 4 mm) coloured green. This output can supply DC voltage or act as a generator of periodic functions.

The outputs S1 and S2 are referenced in relation to zero V accessible on the safety socket (diameter 4 mm.) coloured black located just by the side (mass).


## 4.4 Synchronisation input

This input is accessible on the top of the console, on the safety socket (diameter 4 mm.) coloured red. It is referenced in relation to 0 V accessible on the safety socket (diameter 4 mm.) coloured black located just by the side (mass).

## 4.5 Connection to the computer

Connection to the computer is by the USB socket accessible at the rear of the unit, with the USB lead supplied with the console.


Do not unplug the USB connector from the console when acquisition is in progress.

## 4.6 Connecting a Webcam


On the rear left of the console, there are 3 connectors for image input.


Connect your webcam to the input corresponding to the diagram below.

## 4.7 Connection of a fast camera


- Switch off your computer and unplug the power supply of your console
- In addition to the USB connection, connect your console (fig.1) to the IEEE 1394 port of your computer using the required lead (see table)

IEEE port Of your PC	Type of lead needed	Bandwidth
type a (6 pins) 	a to b	400 Mb / s
type b (9 pins) 	b to b	800 Mb / s

- Connect your fast camera to your console (fig. 2) using the cable supplied with the latter.

Fig.2

- Power-on your PC then power-on your console.

## 4.8 Power supply

To operate correctly, the console must always be connected to the mains using the power supply unit.

**Caution:** Use only the power supply unit supplied with the console.

Do not unplug the power supply of the console when acquisition is in progress.

## 4.9 Updating Ex.A.O. software

In order to use your console in data acquisition mode or in order to obtain more functions linked to images which allow synchronised video and data acquisition, it is necessary to have science workshop software version 2.0 at least.

## 5. Technical characteristics

Designation	Option	Value
ESAO® 4 compatible location	Number of channels	4
ESAO® 3 compatible location	Number of channels	2
Analogue inputs on connector	Number of channels Input range	8 ±5 V
Differential analogue inputs	Number of channels Input range	4 ±5 V
Analogue acquisition	Resolution Standard sampling Rapid sampling (1 channel)	12 bits 2 M Ech/s 10 M Ech/s
Bandwidth	attenuation of 1% -3dB	100 kHz 770 kHz
Protection of analogue inputs	Permanent regimen	±250
Synchronisation	External logic Internal analogue	5 V
Triggering	software Standard Post-synchro (depending on software) Post-synchro (depending on software)) Measuring random events	
S1 analogue output	Output range Continuous mode Function generator mode Resolution	0 à 10 V  12 bits
S2 analogue output	Output range Continuous mode Function generator mode Resolution	±5 V  12 bits
Logic inputs-outputs	Individually programmable Recognition of adapters	32 32
Timer for dating logic and analogue events	Resolution Time base Programmable front	24 bits 50 ns
Frequencemeter/Periodmeter	Resolution Minimum time base Maximum time base	16 bits 100 ns 100 ns
Events meter	Resolution	16 bits
Onboard memory		512 ko

## 6. Recommended configuration

The ESAO® Visio Plus console operates with Microsoft Windows 98 SE and superior.

The recommended minimum hardware configuration is: Pentium III 800 MHz, RAM: 256 Mb

Use of fast cameras requires a much more recent configuration. In this case, you are advised to contact us to check your configuration is suitable.

## 7. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.  
For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
468 rue Jacques Monod  
CS 21900  
27019 EVREUX CEDEX FRANCE

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

**Free service\***

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux