Ex.A.O.

Adaptateur ESAO®

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0 825 563 563*

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1. Description

The sensor + probe CO₂Air/Water set allows measuring the CO₂ concentration in an aqueous or airborne environment. It can be used to carry out a data acquisition.

This sensor can be used with the consoles Visio®, Visio®+, ESAO®4+.

2. First commissioning

- Unscrew the head and then fill it with electrolyte in three quarts
- Tap the head on the work surface to remove the air bubbles
- Make sure that a small quantity of the electrolyte overflows when screwing the head back onto the probe.
- Rinse the probe head with a jet of distilled water
- Wait polarization (paragraph 3)
- Calibrate your probe (paragraph 4)

Important for all products ordered before 11/07/19:

The electrolyte which is included with the sensor + probe CO₂ air/water can only be used for a first use. For a long preservation, we advise you to use the electrolyte ref. 105 175.

3. Polarization

Polarization means that a chemical equilibrium is established in all parts of the glass electrode.

The probe should be polarized in the following cases:

- In initial use
- After changing the electrolyte
- After changing the probe head

The probe is polarized by filling the head with electrolyte and letting it stand with the head pointing downwards for a time period specified in the table below. Polarization does not require electricity.

Case	Polarization time
Initial use	10 hours or about one night
After changing the electrolyte	3 hours
After changing the probe head	3 hours

After the probe is polarized, it must be calibrated.

4. Calibration

The probe matched with the sensor limits the need to calibrate it during each session: it saves a lot of time for the preparation of your practical work.

It is necessary to calibrate your probe in these cases:

- In initial use
- After changing the electrolyte
- After changing the probe head
- For precise measurements

Procedure


- Insert the sensor in the console and connect the console to your computer

- Put the probe in the air
- The sensor requires a 10 minute warm-up


Setting the 0.04%:

- In the software, place the sensor icon on the y axis
- Click on the « Calibration » tab and follow the instructions that appear in the tab
- The red light flashes between red and green. When it is stable on the green, the measurement is stable. Then click on Next.
- A 10 second pre-roll appears and the zero is displayed
- Prepare your calibration solution

Setting the 5.82%:

- Prepare the calibration solution. In a beaker, mix 25 ml of basic calibration solution and 3 ml of acid calibration solution
- ⇒ **Caution, this solution is only reactive during 15 minutes.**
- Install the probe in the solution and wait for the lights stabilization

- Click Next
- After a 10 second pre-roll, your probe is set and you can leave the calibration tab
- If the setting is not correct, you obtain the following message and a help file opens.

You can repeat the procedure after checking the probe.

5. Use

5.1 Connecting the Sensor to the Console

- Place the sensor in front of the console connector
- Slide the sensor following the console guide
- Press until reaching contact with the console
- The sensor requires a 10 minute warm-up time

5.2 In the software

The CO₂ sensor is represented by the icon in the Scientific Workshop general module. In order to carry out an acquisition, slide the icon on one of the axis. For example, for the CO₂ level acquisition as a function of time, slide the CO₂ icon on the y axis and the time icon on the x axis.

6. Storage

Once a year, we recommend you to:

- Unscrew the head of the probe
- Empty the electrolyte from the head
- Rinse the head with distilled water or change it if needed
- Scrub the electrode with 95° ethanol
- Fill the head with electrolyte
- Screw the head of the probe

The probe is stored in open air, with the probe head.

Place the probe with the head facing downwards and protect it from being jarred.

Do not store the probe without electrolyte.

Do not store the probe in water or an aqueous solution.

7. Technical Specifications

Warm-up time	10 minutes
Measurement range in air	0 – 32 % of CO ₂
Measurement range in water	0 – 500 mg/L
Measurement accuracy	1% of the measurement
Response time	90 seconds
Probe service life	>5 years

8. Spare parts

Type	Réf.
Électrolyte for CO ₂ sensor (50 mL)	105 175
Acid solution for CO ₂ probe calibration (50 mL)	107 452
Base solution for CO ₂ probe calibration (250 mL)	107 454
CO ₂ sensor heads	453 152

9. After-Sales Service

The warranty is of 2 years.

For any repairs, kindly contact the **Technical Support on 0 825 563 563**.

The equipment should be returned to our workshops, and for any repairs or spare parts, kindly contact

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux