METHODE ESAO® ADAPTATEUR

ADAPTATEUR pH-METRE ESAO 4

REF. 452 109

I - GENERALITES

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Ordering Code: 452 109

I - GENERAL

This adapter is used to perform pH measurements by means of the ESAO system. It is used with a pH probe (combination electrode P/N 703 170).

The adapter is used to power the probe and to match the measurement signal to the input of the ESAO interface.

Two potentiometers are available for calibrating the pH probe.

II - OPERATION

1. ESAO system

This adapter operated with ESAO 4, i.e.:

- an ESAO 4 interface board
- an ESAO 4 console

Ordering Code: 451 100

Ordering Code: 451 101

2. Adapter connection to ESAO console

The adapter is connected via the connector at the rear of the adapter housing:

- hold the adapter housing horizontally.
- follow the guide slots in the ESAO console.
- insert fully home.

3. pH probe connection

This adapter is used in combination with the following pH probe:

Combination electrode

Ordering Code: 703 170

The probe is connected via the BNC connector on the front side of the adapter.

4. Probe calibration

Two adjustment knobs are accessible to the user on the upper side of the adapter.

The left knob (symbol ) is used to adjust the sensitivity.

The right knob (symbol ) is used for zeroing.

Read the instructions relevant to the sensor.

5. Operation

- ensure that the ESAO board is properly installed in your computer.
- ensure that the experimental console is connected to the ESAO board.
- connect the adapter.
- connect the probe.
- power up the computer.

6. Software

This adapter can be used directly with the dedicated software:

Logi pH

Ordering Code: 000 811

Read the software manual.

To use the adapter with the following general software programs:

Generis

Ordering Code: 000 920

Actilab

Ordering Code: 000 813

Willie

Ordering Code: 000 925

the available characteristics must be as follows:

ESAO code for the adapter:	156
Adapter type:	single adapter
Measuring range:	pH 0 to pH 14

For more details, read the manual relevant to the software used.

III - TECHNICAL DATA

- | | |
|------------------|---|
| - probe input: | BNC connector |
| - response time: | approx. 5 sec |
| - calibration: | one zeroing trimmer
one sensitivity trimmer (slope adjustment) |

IV - MAINTENANCE

For any repair, adjustment or spare parts, contact:

**S.A.V. JEULIN
B.P. 1900
27019 EVREUX CEDEX
FRANCE**