**Méthode ESAO<sup>®</sup>**

***ESAO<sup>®</sup> method***

**Ref :  
452 146**

**Français – p 1**

**English – p 5**

**Version : 6006**

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
## 1 General information

This adaptor helps carry out conductivity measurements with the ESAO4 system.

It is used with the CTM conductimeter probe (part no. 701 329). It is equipped with a manual temperature compensation function.

## 2 USE

### 2.1 ESAO System

This adaptor works with ESAO 4, i.e.

- An ESAO 4 interface card      Part no.: 451 300 or 451 100
  - An ESAO 4 console.      Part no.: 451 101
- Or
- An ESAO 4 Plus console      Part no. : 451 401

### 2.2 Connecting the adaptor to the ESAO console

The adaptor is connected using the connector located at the back of the box:

- Hold the box horizontally
- Follow the guide slots of the ESAO console
- Engage it fully.

## 2.3 Connecting the probe

This adaptor is used with the CTM conductimeter probe (Part no. 701 329), which is connected by two green banana plugs located in front of the adaptor case.

For total protection from external electromagnetic radiation it is recommended that the cord of the probe marked by a white ring be connected to the socket marked by a white circle.

The adaptor supplies power to the probe and adapts the measured signal to the input of the ESAO interface.

## 2.4 Calibrating the probe

*Note on temperature compensation:*

Conductivity is a magnitude whose value depends on the temperature of the solution in which the measurements are made. The influence of temperature on the value measured is of the order of 2% per degree.

Temperature compensation helps bring the conductivity value to the one that will be measured if the solution was at a temperature of 25 °C.

It is thus possible to compare the measurements made at different temperatures.

### *Calibration*

To obtain accurate measurements, it is recommended that the device be calibrated before starting a series of experiments. Calibration is required when the probe is replaced.

Ensure that the temperature compensation adjustment button is set to the temperature of the standard solution adapted to the rated capacity used (refer to the table below).

Immerse the CTM probe in a beaker containing the standard solution.

Press the calibration button until the value displayed is equal to the conductivity value at 25 °C indicated on the flask of the standard solution.

Then adjust the temperature button of the solution studied before making the measurement.

Rated capacity used	Suitable standard solution
0 - 2 mS.cm <sup>-1</sup>	1413 μS.cm <sup>-1</sup> (Part no. : 106 191)
0 -20 mS.cm <sup>-1</sup>	12 880 μS.cm <sup>-1</sup> (Part no.: 106 206)

**Attention:** It is advisable not to change the rated capacity during an experiment. In fact, the device being calibrated for a given rated capacity may display different values from one rated capacity to the other.

## 2.5 Commissioning

- Check that the ESAO card is correctly installed in the computer
- Check that the console is connected to the ESAO card
- Switch on the computer
- Connect the adaptor
- Connect the probe.

## 2.6 Software

This adaptor may be used with the dedicated software:


- Générés conductimetry application
- Generis 4 software
- Generis 5 software (Part no. 000720)

If the adaptor is not automatically recognised then proceed as indicated in 2.7.

## 2.7 Creating the Conductimeter adaptor in the Générés 4 adaptors file.

### 2.7.1 With Generis 4

From the Start menu of Windows, go to Programs and select Générés: the following menu is displayed.

Select the ESAO Adaptor submenu: The following programme is displayed.

In the No. section, enter the code 43 (ESAO4 Conductimeter).

In the name section, enter Conductimeter\_CTM

In the type section, select MONOSENSOR.

The unit is mS/cm.

The rated capacity number is 2.

The rated capacity 1 varies from 0 to 2 mS/cm

The rated capacity 2 varies from 0 to 20 mS/cm

Click on the Save button to store the new adaptor.

### **2.7.2 With Generis 5**

The adaptor is automatically recognised.

## **3 TECHNICAL SPECIFICATIONS**

- Probe input: by two safety sockets  $\varnothing$  4mm
- Rated capacities: 2 mS.cm<sup>-1</sup>, 20 mS.cm<sup>-1</sup>, controlled by the software
- Calibration: in one point by the potentiometer
- Temperature compensation: manual, by a potentiometer from 0 to 50°C  
optimum range from 15 to 35 °C
- Accuracy (excluding the probe): 1% of the entire scale

## **4 After-Sales Service**

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**+33 (0)2 32 29 40 50**

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your  
disposal from Monday  
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries