Ex.A.O.

Data logging

**Réf :
453 064**

Français – p 1

English – p 6

Version : 7112

Bioréacteur 3

Bioreactor 3

1. Généralités

1.1 Principe

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
1 General

1.1 Principle

The bioreactor 3 is a containment designed to study fundamental biological phenomena including respiration, photosynthesis and enzyme activity at cell scale.

It is complete and robust, and its use is justified for Life and Earth Science programs in High School.

All in 1, it combines a **stirrer**, shutters to allow light to enter or to prevent it from entering, a probe support cylinder with captive plugs that can be fitted with three 12 mm diameter probes (O_2 , CO_2 ...) and two 6 mm diameter probes (light, temperature, etc.).


The cell suspension is thus placed in a small defined volume. The stirrer optimises the measurement while keeping the medium homogeneous at all times.

1.2 Detailed description

It comprises:

- a base containing the stirrer ❶
- a metallic lid ❷,
- a vessel with a thermostat-controlled containment ❸,
- 2 metallic shutters ❹,
- a probe support cylinder ❺ with:
 - three 12 mm diameter holes (for the O₂, CO₂, Ethanol, turbidimeter probes),
 - two 6 mm diameter holes (for the thermometer and luxmeter probes)
 - 1 hole for injection
 - 5 captive plugs
- 1 adaptor ring for changing the useful volume (10 or 20 ml) ❻.
- 2 turbulents ❼
- one storage area for shutters ❽
- a potentiometer ❾


Power supply not included. Please command the 12V power supply item no 281 512.

2 Use

2.1 Placement


- 1 – Put the vessel with the experimental medium into place in the base,
- 2 - Put the metallic lid into place.
- 3 – Check that the cylinder ring is in its correct position. You can read the useful volume on the top part of the ring **1**. (The useful volume is 10 or 20 ml, depending on the direction chosen)

- 4 – Put the probes into place in the cylinder, and close the unused holes using the plugs **2**.
- This cylinder is provided with three holes for 12mm diameter probes **3** (O₂, CO₂, ethanol...) and two 6mm diameter holes **4** (thermometer, luxmeter).
- 5 – Put the cylinder into position in the housing of the metallic lid of the bioreactor.

2.2 Hints

The cylinder is provided with a tab **5** at the ring, that is used to put the assembly into place with these probes on the side at the storage bar of the shutters without tilting the probes. Do not forget to connect the stirrer to the mains power supply using the 12V power supply module.

In order to start the motor, it's necessary to :

- put the switch on O position
- connect the bioreactor to the sector
- turn the potentiometer ⑨ as much as possible on the right
- put the switch on I position
- turn gradually the potentiometer ⑨ on the left and put it on the expected speed.


3 Possible manipulations

The main practical work that can be done consists of:

- the study of gas exchanges during photosynthesis
- fermentation of yeasts
- enzymology
- cellular respiration

4 Usage precautions

- Do not clean the metallic part with an abrasive sponge.
- During exposure to light, do not place the lamp at a distance of less than 5 cm for a 40 W bulb or 10 cm for a 100 W bulb.

- Close off unused holes using the plugs integrated into the probe support cylinder.
- Do not hold the lid to carry it.

5 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux