Ex.A.O.[®]

Data logging

**Réf :
453 076**

Français – p 1

English – p 10

Version : 2209

**Sonde CO₂
CO₂ Probe**

1. Introduction

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1. Introduction

The CO₂ probe is a Severinghaus electrode probe.
The electrolyte and the membrane are designed specifically for measuring CO₂. The probe can be used in aqueous or gaseous media.

It will work with the following adapters:

ESAO 3.1 CO ₂ metric adapter	Item no. 452 039
ESAO 4 CO ₂ metric adapter	Item no. 452 121
VTT CO ₂ metric sensor	Item no. 472 102
CO ₂ sensor – advanced	Item no. 482 044

in conjunction with the appropriate interfaces; ESAO 3.1, ESAO 4, VTT/Primo and AirNeXT/Foxy, respectively.

It will also work with a pH meter with an input impedance greater than $10^{12} \Omega$ and equipped with a BNC plug, and which reads the values in mV.

The probe consists of a pH electrode. This electrode is designed especially to measure CO₂ and is immersed in a special electrolyte. The pH of this electrolyte is a function of the CO₂ concentration in the medium in which it is placed. A Teflon membrane selectively permeable to CO₂ separates the electrolyte from the medium to be analyzed.

The probe is connected to the CO₂metric adapters via the shielded cable and the BNC plug.

The following types of experiments can be performed:

- human, animal and plant respiration,
- photosynthesis,
- fermentation.

ATTENTION:


The electrode contains ethylene glycol, which is an irritant. The other chemical compounds are not dangerous.

The pH electrode is made of a special, thin glass and is therefore fragile.


2. How it works

2.1 Description

2.2 Electro-chemical principle of the functioning mechanism

The following reaction takes place in the electrolyte of the CO₂ probe:

The following is the dissociation constant for the reaction:

$$K = \frac{[\text{H}^+][\text{HCO}_3^-]}{[\text{CO}_2]}$$

Because the HCO₃⁻ concentration is high in the electrolyte, it varies little, and one thus obtains [H⁺] = K' [CO₂] = K'' pCO₂ (wherein pCO₂ corresponds to the partial pressure of CO₂). The CO₂ concentration is thus measured by measuring the pH of the electrolyte. The CO₂ probe is not compensated for temperature (in contrast to the oxygen probe). It would therefore be necessary, in theory, to calibrate the probe for each temperature change. In practice, however, using the CO₂ probe at a nearly constant ambient temperature is sufficient.

The CO₂ probe has a response time of ca. 2 min. (10 sec. for the O₂ probe).

Consequently, when you remove the probe from the calibration solution after having adjusted the "slope," it is necessary to wait ca. 5 to 10 minutes for it to display a level lower than 0.2%.

If the probe must be readjusted to zero, wait 15 min. before doing so.

3. Setup

3.1 Assembly

The probe is shipped assembled, but without electrolyte. The probe head is protected by a cap. Before use, simply remove this cap. Be sure to remove it gently or the membrane may be damaged.

If the part clipped on the head removes at the same time as the protection, you have to put it back.

3.2 Filling the probe head

Unscrew the head and then fill it with electrolyte (Item no. 105 175). Tap the head on the work surface to remove the air bubbles. Make sure that a small quantity of the electrolyte overflows when screwing the head back onto the probe. The presence of a large air bubble can be detected by taking a reading with the head pointing upwards and then downwards; the value will fluctuate if there is an air bubble. If this is the case, unscrew the head and top up the electrolyte solution.

When screwing the head back onto the probe, the excess electrolyte will overflow; therefore, rinse the probe head with a jet of distilled water.

3.3 Polarization

The probe is polarized by filling the head with electrolyte (see the preceding section) and letting it stand with the head pointing downwards for a time period specified in the table below.

The probe should be polarized in the following cases:

Case	Polarization time
Initial use	10 hours or about one night
After changing the electrolyte	3 hours
After changing the probe head	3 hours

Note:

Polarization does not require electricity; the operation can be performed without plugging the probe into its adapter or the computer assisted experimental system. For this probe, polarization means that a chemical equilibrium is established in all parts of the glass electrode.

After the probe is polarized, it must be calibrated. See the following section.

4. Calibration

The probe must be calibrated in the following cases:

- Initial use
- After changing the electrolyte
- After changing the probe head
- Whenever a precise measurement is required

Operation Mode:

Notes:

- ▶ Ideally, the calibration should be performed at a temperature approximating that of the experiment planned and in a moderately humid atmosphere (50 to 70% RH). In practice, the ambient conditions of a laboratory or classroom are sufficient for a satisfactory calibration.
- ▶ It should be noted that the level of CO₂ in a classroom of twenty or so students can increase by 0.3% by the end of the period. It is therefore best to perform the calibration in a well-ventilated room.


4.1 With your Primo/VTT or Visio interface

4.1.1 Connect the probe to the adapter and the adapter to the ESAO console.


The adapter must be plugged in 15 minutes before beginning the calibration.

4.1.2 Move the selector on the adapter to the “air” setting.

4.1.3 Turn the slope button to the highest setting (clockwise).

Response of the non-calibrated probe

Response after adjusting the slope button

4.1.4 Using the zero button, set the value displayed to 0.

Press the magnification button and use the zero button  to display 0.0 at +/- 0.1. The value displayed is slightly greater than 0 when you release the magnification button.

4.1.5 Adjust the slope.

Pour 25 ml of calibration solution (sodium hydrocarbonate, item no. 107 454) in a pill bottle or other narrow-necked container, add 3 ml of citric acid (item no. 107 452). Mix and then insert the probe in the solution right away.

While gently agitating, wait about 4 min and adjust the slope button  to obtain the value closest to 5.8% +/- 0.1%.

Note:

After 10 to 15 minutes, this solution is no longer effective. If this time period is exceeded, make up a new solution and adjust the value.

4.1.6 Rinse the probe with distilled water.

4.1.7 Hints and tips

Mark the positions of the buttons on the adapter.

After performing a careful calibration, the positions of the buttons on the adapter can be marked with a felt tip pen. Doing so will enable you make an accurate readjustment quickly without having to recalibrate the probe. These markings are only reusable if the same probe is used with the same adapter.

If you change either the adapter or the probe, the markings are no longer usable.

Calibration frequency

As long as the same adapter is used with the same probe, the calibrated probe will remain correctly calibrated. It is therefore not necessary to recalibrate for each use. The calibration can be checked quickly by measuring the CO₂ in the air exhaled into an erlenmeyer flask or a vessel for measuring human respiration (do not blow directly on the probe head). This reactivity makes it possible to perform qualitative measurements.

4.1.8 With Sérénis software

In order to perform a quick calibration, with the Sérénis software it is possible to calibrate 4 probes at once in the same solution in the Sérénis display mode. After the calibration is completed, attach the adapter-probe assemblies to the other interfaces.

4.2 With your AirNeXT or Foxy interface

- 1- Connect your Foxy interface to the mains.
- 2- Connect the CO₂ probe to the CO₂ sensor.
- 3- Plug in the CO₂ sensor.
- 4- Wait for about 10 minutes until the sensor is ready: during this time a triangle appears above the hourglass and vanish when the sensor is ready.
- 5- Put the probe successively head-up and head-down: the value fluctuates if there is an air bubble. In this case, remove the probe head and fill it up again.
- 6- Select Air by pressing the orange button.
- 7- Press the setup button (violet button) at least for 2 seconds. A triangle blinks above the


pictogram

- 8- When the arrows blink alternately ^①, the probe is ready to be setup.

Probe in the air, set 0 by using the potentiometer ^①, then briefly press the violet button

^② : a triangle blinks above the pictogram  ^③.

9- Making of the calibration solution : pour 20 ml of base solution into a flask.
 Add 3 ml of acid solution.

10- Put the probe into the calibration solution.

The displayed value remains the same.

Caution : the lifetime of the solution is 10 to 15 min. If exceeded, make a new solution.

11- When the arrows blink alternately (0), the probe is ready to be setup.

Turn the potentiometer (1) until the value 5,8 is displayed then press the setup button (2).

12- If the setup is OK a triangle appears above the pictogram « OK ».

In case of probe failure, a triangle appears above the pictogram .

The probe is ready to be used in Air or Water !

13- Rinse the probe with distilled water, dry it (avoid touching the membrane) and let it in the air.

The value displayed must be around 0.

5. Using the probe with a pHmeter

You can use the probe with a high impedance pH meter. Proceed according to the following instructions:

- Set the pH meter to read in mV.
 - Hold the probe in the air vertically with the head pointing downwards. Record the value indicated.
 - Place 25 ml of base solution in a 30 ml bottle and add 3 ml of reactive solution.
 - Shake the bottle to homogenize the mixture.
 - Insert the electrode into the bottle in such a way that there is no air in the bottle (see Figure 1).
 - Record the value indicated by the pH meter (ca. 2 minutes after immersing the probe).
- Remember that the calibration and the measurement must be performed at the same temperature.

It is very important to prepare the calibration solution immediately before beginning the calibration. The solution can only be used once and it is only good for about 10 to 15 min.

Determining the slope:

The calibration solution with 3 ml of acid added corresponds to 89 mg/L of CO₂ dissolved in the water and 5.82% CO₂ in the air. The concentration in the ambient air is 0.03%. The slope S_t of the probe can be determined by the ratio of the difference in voltage read to the difference in concentration:

$$S_t = \Delta U / \Delta (\log p_{CO_2}) \quad \times \quad \text{correction factor}$$

The correction factor = 1 in the air and 1.09 in a solution.

Example: voltage in the air (0.03 % CO₂) = -233 mV
voltage in the solution (5.82 %) = -108 mV

$$S_t = (-233 - (-108)) / (\log 0.03 - \log 5.82) = 54.63 \times 1.09 = 59.5 \text{ mV/decade}$$

The probe signal increases by 59.5 mV when one increases the concentration by a factor of 10.

Measurement:

The concentration of an unknown gas can be calculated by the following formula:

$$\log p_M = \log p_1 - (U_1 - U_M) / S_t$$

wherein: log p₁: logarithm of the partial pressure of CO₂ in the air (log 0.03)

U₁: voltage read in the air

S_t: slope

U_M: voltage read in the unknown gas

log p_M: logarithm of the partial pressure of CO₂ in the unknown gas

example: log p₁ = log 0.03

$$U_1 = -233 \text{ mV}$$


$$U_M = -133 \text{ mV}$$

$$\log p_M = \log (0.03) - ((-233 - (-133))/59.5) = 0.16$$

therefore, the gas contains 10^{0.16} % or 1.44 % CO₂.

6. Malfunctions

If you obtain an unstable reading, make sure there are no air bubbles in the glass electrode. If there are bubbles in the electrode, remove them by shaking the electrode the same way you would a medical thermometer. If the problem persists, change the electrolyte and if necessary the probe head.

If you cannot calibrate the probe (it will not rise to 5.82%, blockage at 10%, etc.), gently clean the tip of the electrode with an absorbant tissue soaked in alcohol, change the electrolyte, and if necessary the probe head.

Note:

The calibration described above is a precise calibration. For qualitative work (variations in CO₂ levels during photosynthesis or respiration), it is not necessary to perform such a precise calibration.

Attention: The electrode must be polarized (see Ch III).

In the event of a persistent blockage of the probe at 10%, remove the probe head and rinse the electrode copiously with distilled water until a response is obtained.

7. Testing the probe

You can test the polarization of the probe by connecting it to a pH meter that reads in mV. The value obtained should be between -180 and -240 mV with the probe in the air.

8. Storage and care

The probe is stored in open air, with the probe
Place the probe with the head facing downwards and protect it from being jarred.

**Do not store the probe without electrolyte.
Do not store the probe in water or an aqueous solution.**

9. Technical specifications

Measurement principle:	<u>Severinghaus</u> principle (electrochemical-potentiometric probe)
Measurement range:	from 0.03% to 100% CO ₂
Operating temperature:	from 5 to 45 °C
Temperature compensation:	NO a variation of 1°C induces a measurement error of 0.03% for 1% and 0.12% for 4%.
Length:	120 mm
Diameter:	12 mm
Construction material of the probe body:	POM

10. Consumable materials

Probe head (with membrane)	Item no. 453 138
Electrolyte (50 mL)	Item no. 105 175
Base solution (250 ml)	Item no. 107 454
Reactive solution (50 ml)	Item no. 107 452

11. After sales service

For any repairs to this apparatus or spare parts during or after the warranty, please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux