**E.x.A.O**

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***Data logging***

*Ecosystem study chamber*

Réf :  
453 158

Français – p 1

English – p 5

Version : 0211

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***Ecosystem study chamber***

## 1. Description

### Principe

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## 1. Description

### Principle

Device for reproducing and studying representative ecosystems on a small scale.

Alone or interconnected, these waterproof transparent chambers are perfect for reproducing a specific environment, desert, forest, meadow, pond etc. Equipped with probe support, they are compatible with ExAO, Redy, Initio and Grove measurement systems for monitoring temperatures and hygrometry, light or oxygen conditions.

Transparent, of limited volume for faster lab work, it is particularly suitable for studying breathing with small animals (mealworms, crickets, stick insects, etc.), plants, mushrooms, etc.

### Composition

- 1 translucent plastic enclosure, int. (L x W x H): 14.2 x 8 x 9 cm
- 1 probe holder
- 1 connector and its O-ring
- 1 rubber stopper
- 1 adapter plug

## 2. Inserting the probe

The sensor is provided with an adapter plug for probe.

### Using the adapter plug

The adapter plug is used to hold any type of probe with a diameter between 12 mm and 22 mm

**The 12 mm diameter probes can be inserted directly.**

For probes with a diameter greater than 12 mm, the cap must be cut to the desired diameter.

**Proceed as follows :** (The example below is provided for a 15.1 mm diameter probe.)

- 1) Insert the probe into the adapter plug

- 2) Locate the stop zone using your fingers

- 3) Cut out, using scissors, the stopper at the mark

- 4) Insert the probe into the cap to check that it is properly maintained

### **3. After-sales service**

The device is under a 2-year guarantee, it must be sent back to our workshops.  
For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
468 rue Jacques Monod  
CS 21900  
27019 EVREUX CEDEX FRANCE

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

**Free service\***

**+33 2 32 29 40 50\*\***

\* Call cost not included.  
\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux