Ex.A.O.

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Datalogging

Primo/VTT sensor

Réf :
472 004

Français – p 1

English – p 3

Version : 8109

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Thermometer sensor -20/120°C

1. Description

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0 825 563 563*

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1. Description

This sensor enables measurement of temperature from -20 to 120°C with the VTT system. The probe used is a 316L stainless steel probe (pharmaceutical quality) of Pt1000 type. This sensor is used with the VTT® and Primo® interfaces.

2. Use

2.1 Connecting the sensor to the console

- Introduce the sensor aligned with the measurement channel chosen.
- Locate the sensor in the console's guide slots.
- Follow the console's guide slot system.
- Locate fully home.


2.2 Use

- Connect the sensor.
- Switch on the console.

2.3 Calibrating the probe

In the event of deviation in measurements, proceed as follows:

- Prepare a recipient, R0, containing a bath of melting ice (at 0°C).
- Prepare another recipient, R100, containing pure, boiling water (at 100 °C)
- Locate the potentiometers 'zeroing access orifices on the front of the sensor (P100 at left and P0 at right).
- Connect the sensor.
- Switch on the console.
- Use the screwdriver supplied with the console to make the following adjustments:
 - 1 - Plunge the probe into R0.
 - Wait for the measurement to stabilize.
 - Set P0 to display 0°C.
 - 2 - Plunge the probe into R100.
 - Wait for the measurement to stabilize.
 - Set P100 to display 100°C.

Repeat points 1 and 2 several times to simultaneously get measurements that are accurate in line with the values stated in paragraph III.

3. Technical features

Length	135 mm
Diameter	6 mm
Measurement range	From -20 to 120°C
Resolution	0,14 °C
Accuracy	±2 °C
Response time in water	About 15 s
Maximum exposure temperature	-50°C/+200°C

4. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux