**Ex.A.O.**

Capteur Primo / VTT

***Data logging***

*Sensor*

Réf :  
472 089

Français – p 1

English – p 5

Version : 8106

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***Air  $CO_2$  Sensor***

## 1. Description

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
**0 825 563 563\***

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## 1. Description

The Air CO<sub>2</sub> sensor measures the CO<sub>2</sub> concentration in air. It can be used for data acquisitions making use of VTT, Primo and Syneo consoles.

## 2. Use

### 2.1 Precautions

- Do not immerse the probe because there is a risk of irreversibly damaging it.
- For better precision, switch the probe power on 5 minutes before starting any manipulation.

### 2.2 Interface

This sensor is used with the VTT, Primo and Syneo consoles.

### 2.3 Inserting the sensor


- Present the sensor facing the connector on the console.
- Slide the sensor following the guide on the console.
- Push until it clicks into place.

## 2.4 Inserting the probe


The sensor is delivered with a probe adapter. It makes it possible to use the probe with the following products:

- Animal breathing chamber (ref: 453 066)
- Human breathing chamber (ref: 453 143)
- Human metabolism chamber (ref: 453 144)

For the animal breathing chamber:

For human breathing chambers:


## 2.5 Adjustment of 0.04 %

- Place the probe in open air.
- Do a long press (a few seconds) on the button positioned on the sensor.
- The LED flashes red for the time of the calibration.
- Once the LED is continuously on in red, calibration is finished.

**Note:** It is possible that the displayed value is slightly different from 0.04%. In this case, start the calibration procedure again.

## 2.6 Data acquisition

The sensor CO<sub>2</sub> is represented by the **CO<sub>2</sub>** icon in the *Atelier Scientifique* general-purpose module. When you want to make an acquisition, slide the icon onto one of the axes. For example, to acquire the CO<sub>2</sub> content as a function of time, slide the CO<sub>2</sub> icon onto the ordinate and the time icon onto the abscissa.

To start acquisition, click on .

## 3. Technical characteristics

<b>Measurement</b>	CO <sub>2</sub> content in air
<b>Measurement range</b>	0-5%
<b>Precision</b>	at 5% of CO <sub>2</sub> : ± 0,03% of CO <sub>2</sub> at 0,2% of CO <sub>2</sub> : ± 0.01% of CO <sub>2</sub> at 0,1% of CO <sub>2</sub> : ± 0,006% of CO <sub>2</sub>
<b>Heating time</b>	30 seconds to 1 minute
<b>Probe useful life</b>	>5 years

## 4. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.  
For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
468 rue Jacques Monod  
CS 21900  
27019 EVREUX CEDEX FRANCE

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

### Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.  
\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux