**Ex.A.O.**

Capteur Foxy/Tooxy/AirNeXT

***Data logging***

*Sensor*

Réf :  
482 034

Français – p 1

English – p 4

Version : 8110

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***Current sensor TRMS***

## 1. Généralités

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**0 825 563 563\***

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## 1. General

This sensor measures effective (TRMS - True Root Mean Square) and instantaneous current intensities up to  $\pm 10$  A, using four calibrations. Measurements are made via a differential input on safety banana jacks: (two red, one black, diameter 4 mm).

## 2. Use

This sensor operates with the Foxy console and can be used:

- in stand-alone mode by displaying the results on the sensor display
- connected to a computer, by displaying the data in the Scientific Workshop software.

## 3. Connection to the Foxy console

The sensor is connected using the connector located on the unit:

- Align with the chosen measuring channel.
- Engage the unit on the console guide.
- Push well home.

## 4. Type of measurement

The sensor can measure the effective (TRMS - True Root Mean Square) value of current intensity, or its instantaneous value.

A button on the sensor unit allows modifying the type of measurement selected, which is then identifiable, in stand-alone mode, by an arrow on the ammeter sensor display.

## 5. Choice of calibration

The four available calibrations are:  $\pm 0.05$  A,  $\pm 0.5$  A,  $\pm 2.5$  A,  $\pm 10$  A.

A button on the sensor unit allows changing the calibration.

For calibrations  $\pm 0.05$  A et  $\pm 0.5$  A, use input mA.

For calibrations  $\pm 2.5$  A and  $\pm 10$  A, use input A.

## 6. Operation in stand-alone mode

### Commissioning:

- Connect the sensor to the console
- Power-on the console


### Display:

An instantaneous or effective (TRMS - True Root Mean Square) value is displayed on the ammeter sensor screen.

### Type of measurement/calibration:

The parameters are modified using two buttons located on the sensor (orange for type, green for calibration).

The selected parameter is indicated by a black arrow: ▲.

## 7. Operation in connected mode

### Commissioning:

- Connect the sensor to the console
- Power-on the console
- Link the console to the computer

### Display:

The ammeter sensor display is inactive, the measurements are displayed directly on the computer screen via the acquisition software.

### Type of measurement/calibration:

These parameters can be modified using the acquisition software and the two sensor buttons are then deactivated.

## 8. Technical characteristics

<b>Measuring range Instantaneous Mode</b>	<ul style="list-style-type: none"> <li>- calibration <math>\pm 0.05</math> A      From <math>-0.05</math> to <math>+0.05</math> A</li> <li>- calibration <math>\pm 0.5</math> A      From <math>-0.5</math> to <math>+0.5</math> A</li> <li>- calibration <math>\pm 2.5</math> A      From <math>-2.5</math> to <math>+2.5</math> A</li> <li>- calibration <math>\pm 10</math> A      From <math>-10</math> to <math>+10</math> A</li> </ul>
<b>TRMS Mode measuring range</b>	<ul style="list-style-type: none"> <li>- calibration <math>0.05</math> A      From <math>0</math> to <math>+0.05</math> A</li> <li>- calibration <math>0.5</math> A      From <math>0</math> to <math>+0.5</math> A</li> <li>- calibration <math>2.5</math> A      From <math>0</math> to <math>+2.5</math> A</li> <li>- calibration <math>10</math> A      From <math>0</math> to <math>+10</math> A</li> </ul>
<b>Accuracy</b>	$\pm 1\%$ of the measurement
<b>Bandwidth</b>	$\geq 20\text{kHz}$
<b>Protection</b>	By a resetting fuse, up to $10$ A

## 9. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.  
For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
468 rue Jacques Monod  
CS 21900  
27019 EVREUX CEDEX FRANCE

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

**Free service\***

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux