Foxy
Capteur

Foxy
Sensor

Réf :
482 038

Français – p 1

English – p 4

Version : 2101

Capteur pression Foxy
Foxy pressure sensor

1 Généralités

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0 825 563 563 *

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1 General

This sensor measures the absolute or relative pressure in relation to atmospheric pressure in non-corrosive gases via a pipe of internal diameter 3 for attachment to an end fitting. This pipe must withstand a pressure of seven bar (crystal type PVC).


2 Use

This sensor can be used with any console in the Foxy range.

2.1 Connection of the sensor to the console

- Position the sensor opposite the connector on the console
- Slide the connector following the console guide
- Push home until clipped in.

2.2 Commissioning

2.3 Use in stand-alone or connected mode

2.3.1 In stand-alone mode

Select the calibration using the calibration button.

The calibration selected is indicated by a triangle in the green strip opposite the calibration.

The absolute pressure is displayed.

If relative pressure calibration is selected, the tare button becomes operational and the this pressure value can be set to zero.

For this, first press down and hold the tare button.

A triangle flashes opposite the symbol .

Wait a few moments then press the tare button again briefly.

The value stabilises around 0.

2.3.2 In connected mode

When connected to a computer, the sensor display goes out. The measured values are visible only on the computer and buttons become inactive.

3 Additional equipment

No additional equipment is required.

4 Technical characteristics

- Calibration-resolution:
 - 7000 hPa 1.75 hPa
 - 800 - 1200 hPa 0.3 hPa
 - ± 200 hPa 0.1 hPa

Remarks on calibration:

Calibration 7000 hPa, given its wide measuring range, is appropriate for the traditional general uses of pressure sensors.

The calibration 800-1200 hPa is specifically provided for measuring atmospheric pressure.

The calibration ± 200 hPa allows making relative measurements and is specifically adapted to TP [Practical Work] where the aim is to follow up pressure variations without taking into account atmospheric pressure, as typically verifying the relation $P=F/S$.

- Accuracy: $\pm 1\% + 2$ LSB of the measurement.

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries