**Ex.A.O.**

Capteur Foxy/Tooxy/AirNeXT

***Data logging***

*Sensor*

Réf :  
482 046

Français – p 1

English – p 5

Version : 8110

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**0 825 563 563\***

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## 1. General

Extremely versatile, the multi-functional Electrophy sensor measures all electrical phenomena in SVT.

It allows measuring both an ECG or the cardiac frequency and also electrical phenomena such as stimulation of a crab nerve or displaying the myotatic reflex.

The Electrophy sensor is supplied with an ECG cable for connection to an individual using self-adhesive electrodes.

## 2. Use

The sensor can be used:

- In stand-alone mode, showing the data on the sensor display.
- Connected to a Foxy console, displaying the data on the computer.

### 2.1 Connection to the console

- Align the unit with the chosen measuring channel.
- Engage the unit on the console guide.
- Push well home.

### 2.2 Use in stand-alone mode

Cardiac frequency:

Only the cardiac frequency mode is available in stand-alone mode.

Connect the cable supplied to the Electrophy sensor. Position three self-adhesive electrodes (ref. 453029) as follows:

- The yellow electrode is located on the left wrist
- The red electrode is located on the right wrist
- The violet electrode is located on the left ankle.

Remark: If measuring with electrodes on the wrist, it is recommended that at least the upper half of the body is at rest, otherwise the system will measure the muscular electrical activity of the arms. For measurement under effort, place the electrodes on the chest.

Positioning on the chest:

Placing the electrodes on the chest produces better quality results. In this case, the electrodes should be positioned as follows:

- The yellow electrode on the left side
- The red electrode on the right side
- The violet electrode is placed on the back, or near the waist line.

After positioning the electrodes, check 'FC' mode is selected, and otherwise, change the mode using the orange button.

For each individual, the gain will be different. If the frequency displays 0, increase the gain until a measurement is obtained and four to five arrows are flashing on the violet band. Adjustment is complete, you can now display the cardiac frequency.

## 2.3 Use in connected mode


In connected mode, the sensor display is inoperative.

### Cardiac frequency

Position the electrode as indicated in paragraph 2.2.

Select the 'cardiac frequency' mode in the 'Calibration' tab.

Adjust the gain using the 'Adjustment' tab. Check the white circle is in the correct zone. Otherwise, adjust it using the cursor located below.

You can then start the readings, generally for a few minutes, to see how cardiac rate increases during effort.

### **ECG**

Position the electrode as indicated in paragraph 2.2.


Select the 'ECG' mode in the 'Calibration' tab.

Adjust the gain in the 'Adjustment' tab.

You can then begin taking readings for ten seconds, for example.

### **Myotatic reflex**

Position the electrodes on the ankle as shown in the diagram below.


Link the ECG cable to the sensor.

Select the 'Electrophy' mode in the 'Calibre' tab.

Connect the 'Martoreflex' to the grounding terminals (black) and the synchro terminals (red) on the front panel of the Foxy console.

Parameter synchronisation on the Martoreflex is indicated below.

You can then begin readings over for 100 ms. Recording will start when the Martoreflex hammer strikes.

## **3. Technical characteristics**

- Measuring range: 0 to 220 cpm (strikes per minute)
- Measuring accuracy:  $\pm 1$  strike per minute of maximum value
- Number of decimal places displayed in stand-alone mode: 1

## **4. After-sales service**

The device is under a 2-year guarantee, it must be sent back to our workshops.  
For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
468 rue Jacques Monod  
CS 21900  
27019 EVREUX CEDEX FRANCE

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

**Free service\***

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux