Foxy
Capteur

Foxy
Sensor

Réf :
482 048

Français – p 1

English – p 5

Version : 2101

Capteur chronociné Foxy
Foxy chronocine-time lapse sensor

1 Généralités

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27019 EVREUX CEDEX France
0 825 563 563*

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1 General

This sensor can detect, measure and control logical events using the dedicated Scientific Workshop module: Chronoméca.

2 Use

This sensor operates with the Foxy console and can be used:

- In stand-alone mode, displaying the results on the sensor display,
- Connected to a computer, visualising data on a screen.

3 Connection to the Foxy console

The sensor is connected using the connector located on the unit:

- Align the unit with the chosen measuring channel
- Engage the unit on the console guide
- Push well home

4 Type of measurement

Equipped with a Chronociné sensor/fork (ref. 453026) or Chrono ESAO (ref/ 453008), it allows:

- Detecting events
- Dating events
- Measuring average or instantaneous speeds

A button on the sensor unit allows, in stand-alone mode, selecting a sensor/fork and displaying the corresponding value.

Supplied with a Y-shaped adaptor, the unit can be combined with 1 or 2 Chronociné sensors (or forks).

5 Operation in stand-alone mode

Commissioning:


- Connect the sensor to the console
- Power-on the console

Display:

After linking the forks to the sensor, the result is read directly on the Chronociné sensor display. Automatic reset is performed at the end of 60 seconds.

To check satisfactory operation of the forks, simply pass an object in front of the infrared beams of the forks and observe appearance of the detection indicators ▼ in the violet zone of the display.

The result displayed on the sensor corresponds to a given fork, indicated by a ▲ in the orange zone.

Selection of sensor/fork:

Observation of the results corresponding to another fork is achieved by pressing the orange button on the sensor.

6 Operation in connected mode

Commissioning:

- Connect the sensor to the Foxy console
- Power-on the console
- Link the console to the computer
- Start the Physics-Chemistry Scientific Workshop then the dedicated Chronoméca module


Display:

The Chronociné sensor display is inactive, the measurements are displayed on the computer screen via the Chronoméca module, after acquisition of an event.

3 forks detected

Positioning of
forks

Selecting a sensor/fork:

In this mode, there is no need to make a selection. All the data for each sensor/fork is taken into account by the dedicated Chronoméca module.

7 Technical characteristics

Calibration	None
Resolution	0.1526 mm/s
Display format	#####

8 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries