Ex.A.O.

Capteur Foxy/Tooxy/AirNeXT

Data logging

Sensor

Réf. :
482049

Français – p 1

English – p 3

Version : 8110

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0 825 563 563*

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1 General

This sensor, with an integral or remote microphone, allows acquiring an audible signal via both microphones and measuring its intensity in decibels using the integral microphone. With a remote microphone, the audible signal is measured exclusively in UA (arbitrary unit).

In this case the function dB(A) is deactivated and the function mode button is inactive.

2 Use

This sensor can be used with any console in the Foxy range.

2.1 Connecting the sensor to the console

- Position the sensor opposite the console connector.
- Slide the connector following the console guide.
- Push home until clipped in.

2.2 Commissioning

Select the micro or sonometer function using the orange button.
Then adjust the level by turning the knob.

2.3 Use in stand-alone or connected mode

2.3.1 In stand-alone mode

The value displayed is that of the sound level in dB.

2.3.2 Connection to a computer

The sensor display goes out. The values measured are visible only on the computer and the buttons become inactive.

The principle of measuring is the same as for stand-alone mode, simply follow the instructions on the computer screen.

3 Additional equipment

No additional equipment is necessary excluding the attaching system.

4 Technical characteristics

Integral electret microphone or external microphone connected by jack, diameter 2.5mm
- sound wave form
- sound intensity from 30 to 120 dB (A)
The default function is: dB(A).

Resolution (12 bits/ 4000 points):

0.0225 dB

0.25 u.a

Accuracy

± (1% +2 LSB) of the value

With the external microphone the function dB (A) is deactivated.

5 After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux