**Foxy**  
Capteur

**Foxy**  
*Sensor*

Réf :  
**482 050**

Français – p 1

English – p 6

Version : 2101

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0 825 563 563 \*

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## 1 General and description


This adaptor allows measuring the turbidity of a suspension and deducing the particle concentration. By extension, it allows calculating the cellular concentration of a solution.


## 2 Operating principle

The intensity of light reduces as it passes from one medium to another medium with a refraction index different from that of the initial medium. This phenomenon is known as dispersion and is observed by the techniques of turbidimetry and nephelometry.

In fact, turbidimetry measures the intensity of the light transmitted, whereas nephelometry measures the intensity of light dispersed in a given direction (perpendicular to the emission of light for the Foxy sensor).

The quantity of light diffused perpendicularly is proportional to the quantity of material in suspension, within the limits of a given concentration (20g/L for yeasts).

### 3 Use

This sensor operates with the Foxy console and can be used:

- In stand-alone mode, showing the results on the sensor display,
- Connected to a computer, displaying the data on the computer.

### 4 Connection to the Foxy console

The sensor is connected using the connector located on the unit:

- Align the unit with the chosen measuring channel
- Engage the unit on the console guide
- Push well home

## 5 Type of measurement

This sensor measures the cloudiness of a liquid from 0 to 100%, and is particularly suitable for monitoring cell growth.

## 6 Operation in stand-alone mode

### Commissioning:


- Connect the sensor to the console
- Power-on the console

### Calibration:

To calibrate the sensor, hold down the zero adjustment button (violet); the arrow ▼ above the "0%" flashes.

Then insert the turbidity meter probe in the standard reference solution (solution 0%, usually distilled water **in a container protected from light**), then briefly press the same violet adjustment button; the arrow above "0%" goes out and the screen displays "0".

The sensor is now calibrated.

### Display:

The turbidity values are displayed directly on the sensor screen.

## 7 Operation in connected mode

### Commissioning:

- Connect the sensor to the console
- Power-on the console
- Link the console to the computer

### Calibration:

*Remark: if the sensor has already been calibrated in stand-alone mode, it is not necessary to recalibrate using the software.*

The zero adjustment button is inactive when the console is connected to the computer.

An “Adjustment” tab in the acquisition software allows calibrating the sensor using a zero solution (distilled water in a **container protected from light**).

Before clicking on adjustment:

- Immerse the probe in distilled water in a recipient protected from light (figure 1)
- Enter the calibration value using the keyboard, here 0 % (figure 2)
- Confirm calibration by clicking the “adjust” button (figure 3)

Figure 1


Figure 2



Figure 3

### Display:

The turbidity meter sensor display is inactive, the measurements are displayed directly on the computer screen via the acquisition software.

## 8 Technical characteristics

<b>Measuring range</b>	0 to 100 %
<b>Resolution</b>	0.025 %
<b>Accuracy</b>	$\pm (1 \% + 2 \text{ LSB})$ of the value
<b>Calibration</b>	At 1 point with the standard reference solution

## 9 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
Rue Jacques Monod  
BP 1900  
27019 EVREUX CEDEX FRANCE

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your  
disposal from Monday  
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries