Foxy
Capteur

Foxy
Sensor

Réf :
482 055

Français – p 1

English – p 4

Version : 2101

Capteur Sismo Foxy

Foxy Sismo sensor

1 Généralités

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0 825 563 563 *

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1 General

The Sismo sensor records the wave generated by a shock or a mechanical breakage.

The signal is measured in an arbitrary unit, the SWA or Seismic Wave Amplitude.

2 Use


The Sismo sensor is used only in connected mode, that is by displaying data on a computer via the Foxy console.

2.1 Connection to the console


- Align the unit with the chosen measuring channel.
- Engage the unit on the console guide.
- Push well home.

2.2 Connected mode

Connect the Foxy console to the computer using a USB connection.
The software starts automatically.

Slide the Sismo icon on one of the points on the ordinate axis.
In the parameters for the Sismo sensor, an “Adjustment” tab allows increasing or reducing the signal. For example, if on seeing the results, the signal is peak-shaved, reduce the gain.
The gain can also be modified using the adjustment button on the Sismo Adaptor.

3 Examples of operations

3.1 Comparison of propagation of a wave according to the type of support medium

With 2 Sismo sensors, it is possible to calculate the speeds. On the surface under investigation position the first red sensor, used as a trigger, and then the second, to receive the signal.

Place the sensors on different surfaces or different support media, such as the seismology materials (concrete Ref. 189083 or wood Ref. 507051).

You can also compare the behaviour of the wave on different support media with a single sensor, by synchronising acquisition using a Martoreflex for example (Ref. 453068).

3.2 Influence of temperature on the seismic wave

Thanks to the seismic and temperature scale model (ref. 507043), with a single sensor it is possible to determine whether the temperature influences propagation of the seismic wave.

By adding a temperature sensor (ref. 482101), at the same time, you can demonstrate the variation in speed of the seismic wave, versus the temperature of the support medium.

4 Technical characteristics

- Piezo-electric cell
- Calibration: ± 100 Seismic Wave Amplitude
- Adjustment of the gain by code wheel
- Use in connected mode only

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries