Ex. A. O.

Capteur

Data logging

Sensor

Réf :
482 070

Français – p 1

English – p 4

Version : 6111

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
0 825 563 563*

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1. Description

The Specific Ion/pH sensor is made to measure:

- the pH of a solution using a pH electrode (not provided)
- the temperature of the solution (supplied)
- the concentration via a specific ion electrode (not provided)

2. Use

2.1 Datalogging

This sensor is used with Foxy®, AirNeXT® and Tooxy® interfaces.

2.2 Inserting the sensor

- Place the sensor facing the interface connector.
- Slide the sensor following the interface guide.
- Push until it clicks.

2.3 Inserting the electrode

The pH or specific ion electrode is fitted to the BNC input.
The temperature probe is inserted into the jack socket.

2.4 Use in measurement mode

Once the sensor is inserted into the interface, the value is automatically displayed on the sensor.

Measurement change:

To move from one measurement to another, press the orange button.

Calibration:

To calibrate the electrode:

- Press and hold the purple calibration button, a flashing arrow appears above the note "solution 1".
- Place the electrode in the first calibration solution.
- Enter the solution value using the knurled knob.
- Press the calibration button briefly, a flashing arrow appears above the note "solution 2".
- Place the electrode in the second calibration solution.
- Enter the calibration solution value using the knurled knob.
- Press the button briefly to end calibration.

The temperature probe allows manual compensation of the value during calibration. The temperature measurement must be taken in advance.

2.5 Use with an interface


The buttons are inactive when the sensor is connected to an interface. Adjustments are made via the software.

Two icons appear on the scientific workshop software interface:

-  to measure pH and ions
-  to measure temperature

To take a measurement, drag an icon onto one of the axes depending on the desired acquisition. For example, to measure pH based on a value entered on the keyboard, place

the  icon on the X axis and the  icon on the Y axis:

Measurement change:

To move from pH measurement to ion measurement:

- Place the  icon on the X axis.
- In the adjustment tab, click on "pH" or "Voltage".

Calibration:

Once the measurement type is validated:

- Click on the "Adjustment" tab.
- Click on "Start".
- Place the electrode in the first calibration solution.
- Enter the solution value in the box.
- Click on "Adjust".
- Place the electrode in the second calibration solution.
- Enter the calibration solution value in the box.
- Click on "Adjust".

The system is calibrated.

The temperature probe allows manual compensation of the value during calibration. The temperature measurement must be taken in advance.

3. Technical specifications

pH measurement

Measurement range	0 to 14 u.pH
Accuracy	0.1 u.pH

Voltage measurement

Measurement range	-2000 mV to 2000 mV
Accuracy	1% of measurement

Temperature measurement

Measurement range	20°C/120°C
Accuracy	0.5°C

4. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

+33 (0)2 32 29 40 50