Ex.A.O.

Capteur Foxy / Tooxy / AirNeXT

Data logging

Sensor

Réf :
482 076

Français – p 1

English – p 5

Version : 0211

Capteur CO₂ Air
Air CO₂ Sensor

1. Description

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
0 825 563 563*

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1. Description

The Air CO₂ sensor measures CO₂ concentration in the air. It can be used for data capture or a direct reading via its display.

2. Use

2.1 Warning

Do not immerse the probe in liquid or it may be damaged permanently.

For a better precision, switch on the probe 5 minutes before any manipulation.

2.2 Interface

This sensor is used with Foxy®, AirNeXT® and Tooxy® consoles.

2.3 Inserting the sensor

- Place the sensor facing the console connector
- Slide the sensor following the console guide
- Push until it clicks

2.4 Inserting the probe

The sensor is provided with an adapter plug for a probe.

Using the adapter plug

The adapter cap is used to hold any type of probe with a diameter between 12 mm and 22 mm

The 12 mm diameter probes can be inserted directly.

For probes with a diameter greater than 12 mm, the cap must be cut to the desired diameter.

Proceed as follows : (The example below is provided for a 15.1 mm diameter probe.)

- 1) Insert the probe into the adapter plug

- 2) Locate the stop zone using your fingers

- 3) Cut out, using scissors, the stopper at the mark

- 4) Insert the probe into the cap to check that it is properly maintained

It allows the probe to be used with the following products

- Animal respiration device (ref: 453 066)
- Human respiration device (ref: 453 143)

2.5 Use in measurement mode

- The sensor has a 1 minute warm-up period.
- During this time, the display is inactive. Once the warm-up period is complete the CO₂ rate appears on the screen.

0.04% adjustment:


- Push and hold the calibration button
- Wait until the value appears on the screen
- Push the button briefly to confirm the calibration

Note: it's possible that the posted value differs slightly from 0.04 %. In this case, begin again the procedure of calibration.

2.6 Use with a console

The buttons are inactive when the sensor is connected to a console. Adjustments are made via the software.

The CO₂ sensor is represented by the **CO₂** icon in the Data Science Studio general module. To complete a capture, drag the icon along one of the axes. For example, to capture the CO₂ rate based on time, drag the CO₂ icon along the Y axis and the time icon along the X axis.

0.04% adjustment:

- Place the sensor icon on the Y axis
- Click on the "Adjustment" tab
- Click on Start
- Enter 0.04% in the box
- Wait until the red rectangle turns green.
- Click on "Adjust"

3. Technical characteristics

Warm-up period	30s - 1 minute
Measurement range	0 - 10% of CO ₂
Measurement accuracy	at 10% of CO ₂ : $\pm 0,8\%$ of CO ₂ at 5% of CO ₂ : $\pm 0,03\%$ of CO ₂ at 0,2% of CO ₂ : $\pm 0,01\%$ of CO ₂ at 0,1% of CO ₂ : $\pm 0,006\%$ of CO ₂
Probe useful life	>5 years

4. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux