Ex. A. O.

Capteur

Data logging

Sensor

Réf :
482 078

Français – p 1

English – p 5

Version : 6111

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Heart rate Sensor

1. Description

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0 825 563 563*

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
1. Description

The heart rate sensor performs two functions:

- It measures the heart rate in beats per minute (bpm)
- It provides information about the shape of the cardiac signal

The measurement is made using the diode-photodiode pair on each side of the thimble. The thimble is fixed to the end of the index finger.

2. Use

2.1 Precautions

This sensor is not designed for a use during prolonged exercise. For this type of use, we recommend cardio frequency meters to be fixed to the wrist or chest.

2.2 Interface

This sensor is used with the Tooxy®, Foxy® and AirNext® interfaces.

2.3 Inserting the sensor

- Connect the thimble to the sensor using the Jack connector.
- Present the sensor facing the console connector.
- Slide the sensor following the guide on the console.
- Push until it clicks into place.

2.4 Use in measurement mode

- Place the thimble on the index finger without forcing.
- A value is displayed on the screen after a few seconds.
- The intensity of the measured signal is indicated by a series of triangles (between 1 and 4), at the bottom right of the display.
- Adjust the position of the thimble on the finger so as to have at least two triangles.

2.5 Use with a console


Two channels are available during a use with a console.

2.5.1 Measuring the heart rate

- Place the "Heart rate" icon  onto the ordinate.
- Place the time icon onto the abscissa.
- Set the time parameter as required.
- Start the acquisition by clicking on the green light.

Note: A 1 point per second sampling period is sufficient for the heart rate.

2.5.2 Cardiac signal shape

- Place the "Pulse" icon  on the ordinate.
- Place the time icon on the abscissa.
- Set the time parameter as required.
- Start the acquisition by clicking on the green light.

Note: For the shape of the cardiac signal, it is recommended that points should not be separated by more than 50 ms.

3. Example experiments

3.1 Heart rate variation during exercise

The impact of an effort on the heart rate can be observed during a short and intense activity.

- Place the thimble on the patient's index finger.
- Set acquisition parameters by placing the "Heart rate" icon on the ordinate and the time on the abscissa.
- Set time parameters as follows:
 - o Duration 3 minutes
 - o Number of points 181 points

The experiment takes place in 3 steps:

- 1 rest phase for 60 seconds
- 1 exercise phase for 10 seconds during which the patient performs a series of flexing movements at a sustained rate
- 1 rest phase during the remainder of the acquisition


It is recommended that the patient should minimise movements of the hand on which the thimble is fixed during the exercise phase.

3.2 Determination of the average heart rate from the signal

The periodic nature of the cardiac activity can be demonstrated from the shape of the cardiac signal, the period of the signal can be measured and the average heart rate can be determined.

- Place the thimble on the patient's index finger.
- Set acquisition parameters by placing the "Pulse" icon onto the ordinate and the time icon onto the abscissa.
- Set time parameters as follows:
 - o Duration 15 seconds
 - o Number of points 301 points
- Start the acquisition.

A curve like that presented below is obtained:

The average period is measured by measuring the time for ten periods and dividing by ten. The frequency is determined by taking the inverse of the period.

4. Technical characteristics

Heart rate

Range	from 0 to 220 bpm
-------	-------------------

Cardiac signal shape

Range	± 100 u. a.
-------	-----------------

Technology

Optical system

5. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux