Ex.A.O

Capteur Foxy / Tooxy / AirNeXT

Data logging

Sensor Foxy / Tooxy / AirNeXT

**Réf :
482 082**

Français – p 1

English – p 8

Version : 5207

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09 69 32 02 10 - Appel non surtaxé

You have just purchased an oxygen probe and we thank you for your confidence. This document aims at guiding you in using your probe.

1. Description

This system consists of four elements:

- The sensor
- The connecting cable
- The probe body
- The probe head

This probe is ready for use and operates with the Tooxy®, Foxy® and AirNeXT® interfaces and take measurements in air and water.

2. Use Precautions

Do not touch, scratch or damage the silicon part of the probe head at the risk of having to change it.

To avoid any damages, this part is protected by range-gating.

3. First Use

Your probe comes ready to start.

- Remove the probe protective cap.
- Insert the probe into a channel of an interface
- Wait a few seconds for the measurement to stabilize
- If the value that appears on the sensor screen is different than 20.9 %, then manually set the slope, as follows:
 - Tap and hold the purple button (3 seconds)
 - Once 20.9 % appears, press the purple button briefly
 - The slope of your probe is set.

Your probe is ready to operate

Note: this setting can also be done in the software.

4. Measure

4.1 Using the Probe in Standalone Mode:


The probe takes measurements in air and water. To go from one mode to another:

- In standalone mode, press on the orange button
- In the software, select « Water » in the sensor setting


In the water, in case of outliers, an air bubble may be present under the probe head. To remove it, tap on the probe body.

4.2 Using the Probe with the Scientific Workshop Software:


- Run the Standard Scientific Workshop software


- Drag and drop the oximeter sensor in ordinate:


- The tabs specific to the sensor appear:

- In the caliber tab, select Air:


- Then click on the Measure tab:


- If the value displayed is very different from 20.9 %, click on the Regulation tab

4.2.1 Slope Calibration with Computer

To adjust the slope, you only need to follow the instructions of the Scientific Workshop software:


After fully acquainting yourself with the information, click on Start for the following window to appear:


The screenshot shows the 'Oxymeter' software window. It has a tabbed interface with 'Caliber', 'Quantity', 'Measure', 'Color', and 'Regulation' tabs. The 'Caliber' tab is active. The text inside says 'Enter the value of calibration, then click 'Adjust''. Below this, there is a text input field containing '20.9' and a label '%O2'. At the bottom, the 'Adjust' button is highlighted with a red rectangular border. A mouse cursor is visible near the bottom right corner of the window.

The « Adjust » button is surrounded with a red edging which indicates that the system is verifying the data provided by the probe.
Let this step happen, until a green edging appears:

This screenshot is identical to the previous one, but the 'Adjust' button is now highlighted with a green rectangular border, indicating that the system has accepted the calibration value.


You can then click on the « Adjust » button

The screenshot shows the 'Oxymeter' window after clicking 'Adjust'. The text inside now says 'Your probe is calibrated!'. The 'Adjust' button is no longer visible, and the window has a clean, confirmed state.

The confirmation window for the slope adjustment appears.

You can control the setting by clicking on the « Measure » tab to notice that the value measured by the probe corresponds to the O₂ rate in the Air, being 20.9 %.

You may now use your probe.

4.2.2 Slope Calibration in Standalone Mode

The slope calibration can also be made directly in standalone mode, thanks to the sensor screen.
See the detailed procedure in chapter 6.3 below.

5. Maintenance - Storage

Maintenance: After each use, rinse the probe head with distilled water.

Storage:

For short storage periods (e.g. between 2 PT sessions): Replace the protective cap. Be sure to keep the foam in the cap moist.

For long storage (during summer vacations): Replace the protective cap. Keep the foam in the cap dry.

We advise you to moisten the foam the day before the next TP.

6. Changing the Probe Head

The probe head is designed to function 2 years. Beyond this period, it is preferable to replace it with a new one (ref 453153).

6.1 Replacement

- Unscrew the existing probe head
- Screw the new probe head

Once this operation is finished, a full calibration (zero adjustment + slope) of the probe is necessary:

6.2 Zero Calibration

- Insert the sensor into a channel of an interface
- On the sensor, go to « Eau » mode (meaning “water” mode) by clicking on the orange button
- Immerse the probe into a zero solution (saturated solution of anhydrous sodium sulfite) during 2 to 5 minutes
- Tap and hold the purple button for 10 seconds
- Wait 30 seconds
- Click briefly
- Remove the probe from the solution and rinse it
- The zero adjustment is finished.

6.3 Slope Calibration

- Insert the sensor into a channel of an interface.
- Leave the probe in the open air
- On the sensor, go to « Air » mode by clicking on the orange button
- On the sensor, tap and hold the purple button
- Wait 10 seconds
- Click briefly on the sensor
- The slope adjustment is finished.

7. Accessories

O₂ probe head for optical probe ref. 453 153.

Changing it may only take place in case of repair, if the probe fails to respond or the head is broken. (See chapter 6)

This head does not contain a lining or an electrolyte (optical functioning and not chemical).

8. After-Sale Services

The warranty is of 2 years.

For all settings, contact our **Technical Support** on **+33 2 32 29 40 23**.

The material must be returned to our workshops and for all repairs or spare parts, kindly contact:

JEULIN – S.A.V.
468 Jacques Monod Street
CS 21900
27019 EVREUX CEDEX France

+33 2 32 29 40 23