Ex.A.O.

Capteur Foxy / Tooxy / AirNeXT

Data logging

Sensor

Réf :
482 085

Français – p 1

English – p 6

Version : 2209

Capteur + sonde CO₂ Air/Eau
CO₂ Air/Water Sensor + Probe

1. Description

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0 825 563 563*

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1. Description

The sensor + probe CO₂ Air/Water set allows measuring the CO₂ concentration in an aqueous or air environment. It can be used to carry out data acquisition, or in direct reading thanks to its display.

This sensor is used with the Foxy®, AirNeXT® and Tooxy® interfaces.

2. First commissioning

- Unscrew the head and then fill it with electrolyte in three quarts
- Tap the head on the work surface to remove the air bubbles
- Make sure that a small quantity of the electrolyte overflows when screwing the head back onto the probe.
- Rinse the probe head with a jet of distilled water
- Wait polarization (paragraph 3)
- Calibrate your probe (paragraph 4)

Important for all products ordered before 11/07/19:

The electrolyte which is included with the sensor + probe CO₂ air/water can only be used for a first use. For a long preservation, we advise you to use the electrolyte ref. 105 175.

3. Polarization

Polarization means that a chemical equilibrium is established in all parts of the glass electrode.

The probe should be polarized in the following cases:

- In initial use
- After changing the electrolyte
- After changing the probe head

The probe is polarized by filling the head with electrolyte and letting it stand with the head pointing downwards for a time period specified in the table below. Polarization does not require electricity.

Case	Polarization time
Initial use	10 hours or about one night
After changing the electrolyte	3 hours
After changing the probe head	3 hours

After the probe is polarized, it must be calibrated.

4. Calibration

The probe matched with the sensor limits the need to calibrate it during each session: it saves a lot of time for the preparation of your practical work.

It is necessary to calibrate your probe in these cases:

- In initial use
- After changing the electrolyte
- After changing the probe head
- For precise measurements

Your probe can be calibrated with two methods: on the sensor or with the software.

Preparing the calibration solution:

In a beaker, mix:

- 25 ml of basic calibration solution
 - 3 ml of acid calibration solution
- ⇒ **Caution, this solution is only reactive during 15 minutes.**

4.1 Use in Stand-alone Mode


- Place the sensor in the interface (do not connect it to your computer)
- Place the probe in the air
- Press and hold the calibration button on "air"
- Wait until the value appears on the screen
- Turn the knob to display « 0.0 »
- Short press to validate the calibration of this first point = the sensor then automatically goes to the next setting
- Place the probe in the 5.8% solution
- Turn the knob to modify the value
- Short press again on the button to validate the calibration.
- Rinse the probe with distilled water, dry it (avoid touching the membrane) and let it in the air. The value displayed must be around 0.

4.2 Using with an interface

- Place the sensor in the interface and connect it to your computer
- In the software, place the sensor icon on the y axis

Setting the 0.04%:


- Place the probe in the air
- Click on the « Setting » tab
- Select 0.04 % in the box

- Wait about 1 minute
- Click on « Set »

Setting the 5.82%:

- Place the probe in the 5.8% solution
- Select 5.82 in the box

- Wait about 1 minute
- Click on “set”
- Rinse the probe with distilled water, dry it (avoid touching the membrane) and let it in the air. The value displayed must be around 0.

5. Use

5.1 Connecting the Sensor to the Interface

- Place the sensor in front of the interface connector
- Slide the sensor following the interface guide.
- Press until it clicks.
- The sensor requires a 10 minute warm-up time.


5.2 Use in Stand-alone Mode

During the 10 minute warm-up time, the display is inactive. Once the warm-up time finished, the CO_2 level appears on the screen.

5.3 Using with an interface

The buttons are inactive when the sensor is connected to an interface. Settings are thus made via software:

The CO_2 sensor is represented by the icon CO_2 in the Science Studio general module. In order to carry out an acquisition, slide the icon on one of the axis. For example, for the CO_2 level acquisition as a function of time, drag and drop the CO_2 icon on the y axis and the time icon on the x axis.

6. Technical Specifications

Warm-up time	10 minutes
Measurement range in air	0 – 32 % of CO_2
Measurement range in water	0 – 500 mg/L
Measurement accuracy	1% of the measurement
Response time	90 seconds
Service life of the probe	>5 years

7. Storage

Once a year, we recommend you to:

- Unscrew the head of the probe
- Empty the electrolyte from the head
- Rinse the head with distilled water or change it if needed
- Scrub the electrode with 95° ethanol
- Fill the head with electrolyte
- Screw the head of the probe

The probe is stored in open air, with the probe head.

Place the probe with the head facing downwards and protect it from being jarred.

Do not store the probe without electrolyte.

Do not store the probe in water or an aqueous solution.

8. Spare parts

Type	Réf.
Électrolyte for CO_2 sensor (50 mL)	105 175
Acid solution for CO_2 probe calibration (50 mL)	107 452
Base solution for CO_2 probe calibration (250 mL)	107 454
CO_2 sensor heads	453 152

9. After-Sales Service

The warranty is 2 years.

For any repairs, kindly contact the **Technical Support on +33 (0)2 32 29 40 50.**

The equipment should be returned to our workshops, and for any repairs or spare parts, kindly contact

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux