# ExAO

Capteur Foxy / Tooxy / AirNeXT

## ***Data logging***

Sensor

Réf :  
**482 094**

Français – p 1

English – p 5

Version : 3212

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**0 825 563 563\***

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- General

This sensor allows you to directly digitize angle, speed, and distance measurements for your mechanical experiments.


This sensor is specifically adapted for the following experiments:

- Simple pendulum
- Vertical oscillation
- Study of circular movements
- Displacement measurements
- 

The pulley allows you to directly attach accessories such as the weighing pendulum.

A square section extension, suitable for the 6 mm diameter axle, is provided to allow the pulley to be moved from a support leg if necessary. (Example: use with the simple and heavy pendulum kit).

## • Use

This sensor has to be used with any console in the Foxy range.

### 1.4 Connection to the Foxy console

- Place the sensor opposite the console connector.
- Slide the connector following the console guide.
- Push home until it clips in.

### 1.5 Commissioning

Before starting the measurements, it is important to perform the mechanical zero adjustment. You can do this directly on the pulley or on the sensor in autonomous or connected mode.

Place the pulley on its support.

Attach any accessories, wait until there are no more movements.

Loosen the fastening slightly and orient the pulley so that the index is as close as possible to the zero of the screen printing.

and the pulley is no longer moving, stay.


calibration  
button


## 1.6 Use in stand-alone or connected mode

### 1.6.1 In stand-alone mode

#### Adjustment:

Select the calibration and mode then adjust zero:

A long press on the adjustment button: Segment  flashes.

A short press on the adjustment button: Taring OK, segment  goes out.

After adjusting zero, the value displayed is the instantaneous measured value.

### 1.6.2 When connecting to the computer

The sensor display goes out. The values measured are visible only on the computer screen and the buttons are inactive.

The measuring principle is identical as for stand-alone mode, simply follow the instructions on the computer screen.

- **Technical characteristics**

- 3 modes of measurement :
  - o Angle (°)
  - o Distance (m)
  - o Angular velocity (rpm or rad/s)

- Angle mode:

**Caliber:  $\pm 180^\circ$  /  $\pm 360^\circ$**

Graduations at degree intervals up to  $\pm 45^\circ$

- Velocity mode:

**$\pm 3000$  rpm /  $\pm 314$  rad/s**

- Distance mode:

**Radii: 10mm / 20 mm / 30 mm**

- **After-Sales Service**

The device is under a 2-year guarantee, it must be sent back to our workshops.  
For any repairs, adjustments or spare parts please contact:

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

### Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux