**Ex.A.O.**

Capteur Tooxy/Foxy/AirNeXT

***Data logging***

*Sensor*

Réf :  
482 102

Français – p 1

English – p 4

Version : 8110

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***Luxmeter sensor***

## 1 Généralités

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**0 825 563 563\***

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## 1 General

This sensor is used to measure indoor or outdoor illumination.

## 2 Use


This sensor operates with the interfaces Tooxy, Foxy and AirNeXT in standalone mode (without computer) or in connected mode (with computer).

### 2.1 Connection to the interface

- Bring the sensor in front of the console connector
- Slide the connector along the console guide rail
- Insert until it clicks into place.

### 2.2 Start up

Choose the measurement range by pressing on the button.

*No calibration is necessary*

## 2.3 Use in standalone mode and in connected mode

In standalone mode, change the range by pressing on the sensor range selection button. The display tab then indicates the selected range.

When the sensor is connected to the computer, the sensor display screen goes off. Measured values are only visible on the computer. The range is chosen by the software.

## 3 Additional equipment

The sensor is supplied with its probe.

## 4 Technical characteristics

Ranges	<ul style="list-style-type: none"><li>○ 0 to 2 kLux</li><li>○ 0 to 200 kLux</li></ul>
Resolution	<ul style="list-style-type: none"><li>○ Range 2 kLux: 0.5 Lux</li><li>○ Range 200 kLux: 50 Lux</li></ul>
Precision	$\pm (1 \% + 2 \text{ LSB})$ of the value
Measurement time with onboard software	1 s

## 5 After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops. For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
468 rue Jacques Monod  
CS 21900  
27019 EVREUX CEDEX FRANCE

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

### Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.  
\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux