Ex.A.O.

Capteur Tooxy/Foxy/AirNeXT

Data logging

Sensor

Réf :
482 106

Français – p 1

English – p 5

Version : 8110

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0 825 563 563*

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1 General

The Air/Water oxygen meter sensor associated with the O₂ probe ref. 453 052 (not supplied) measures the variation of the O₂ content in air or in water.

The probe is connected to the BNC connector on the sensor.

2 Use

The oxygen meter sensor can be operated:

- In standalone mode, displaying data on the sensor display screen.
- In connected mode connected to a Tooxy console and displaying data on a computer.

2.1 Connection to the console

- Bring the sensor in line with the selected measurement channel.
- Insert the sensor along the console guide rail.
- Push all the way in.

2.2 Connection of the probe

The O₂ probe is connected to the BNC connector on the top of the sensor.

2.3 Calibration

- **Measurement in air:**

By default, the selected measurement medium is air. An arrow is displayed in front of the "Air" and "25 % O₂" readings on the sensor.

This sensor operates with a console, in standalone mode (without computer) or in connected mode (with computer) after being connected to the console.

If you are making a measurement in Air, simply fix the value to 20.9 % O₂. If the displayed value is too far from 20.9 %, make a long press on the calibration button. A flashing arrow appears above the 20.9% reading.

Then adjust the value to 20.9 % using the central potentiometer. Once the value has been adjusted, press the calibration button again. The OK probe segment lights up to indicate that the probe is working normally.

- **Measurement in water:**

If you want to make a measurement in water, press the orange button to change the range. The arrows appear in front of the "Water" and "20 mg/l O₂" readings. The procedure for adjusting the displayed value is exactly the same as for the adjustment for a measurement in air.

Note:

If you want to calibrate the probe as a function of the water temperature, please refer to the instructions for probe ref. 453 052.

3 Technical characteristics

Air mode:

Measurement range: 0 to 25 % O₂

Sensor precision: $\pm (1 \% + 2 \text{ LSB})$ of the value. The final precision depends on the probe and calibration.

In standalone mode, the sensor displays the value of the O₂ content with a 1:10 display precision

Water mode:

Measurement range: 0 to 20 mg/l O₂

Precision: $\pm (1 \% + 2 \text{ LSB})$ of the maximum value at full scale

In standalone mode, the sensor displays the value of the O₂ content with a display precision of 1:10

4 After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux