Ex.A.O.

Capteur Foxy / Tooxy / AirNeXT

Data logging

Sensor

Réf :
482 109

Français – p 1

English – p 4

Version : 8106

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0 825 563 563*

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1 General

The Ventilation sensor is used to study human breathing. This sensor makes 3 measurements, namely the breathing rate, volume and flow rate.

It is composed of the adaptor itself that is fixed onto the Tooxy, Foxy or AirNeXT console and the ventilation sensor.

2 Use

2.1 Connection to the console

- Bring the sensor in line with the selected measurement channel.
- Insert the sensor along the console guide rail.
- Push all the way in.

The Ventilation sensor may be used:

- in standalone mode, displaying data on the sensor display screen.
- in connected mode, connected to a Tooxy console and displaying data on a computer.

The button on the front of the sensor is used to modify the sensor function.

There are 3 possible functions:

- The breathing rate in cycle/minute
- The breathing volume in litre
- The breathing rate in litre/s

IMPORTANT

At every measure, put the sensor in the position of use and press the Reset button. The button will blink. Do not move the ventilation adaptor while the measurement is being made, otherwise the values may be disturbed.

2.2 Use and adjustment

In standalone mode, the value is displayed on the adaptor. Choose the required function, and if the displayed value is not equal to zero, press the Reset button on the Ventilation sensor.

After 3 flashes of the red LED, the displayed value returns to zero and you can make your measurement.

In connected mode, the Reset function is done automatically when acquisition starts.

Note:

The sensor display screen is inactive during use in connected mode.

3 Technical characteristics

Breathing rate	
Measurement range	0 to 60 cycles per minute (cpm)
Resolution	0.1 cpm
Measurement precision	± 2 % of the maximum value at full scale
Number of decimals displayed in standalone mode	1
Volume	
Measurement range	± 10 litres and ± 100 litres
Resolution	0.3 mL
Measurement precision	5% of the maximum value at full scale
Number of decimals displayed in standalone mode	2 (± 10 litres) and 1 (± 100 litres)
Breathing flow rate	
Measurement range	±10 L/s
Resolution	0.2 mL
Measurement precision	5 % of the maximum value at full scale
Number of decimals displayed in standalone mode	2

4 After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux