**Ex.A.O.**

Capteur Foxy / Tooxy / AirNeXT

***Data logging***

*Sensor*

Réf :  
482 111

Français – p 1

English – p 4

Version : 2312

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***Particle counter***

## 1. Généralités

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## 5. General

This sensor makes it possible to determine the concentration of fine dust particles in the air. Directive 2008/50/EC defines the tolerance thresholds for the level of PM2.5 and PM10 per m<sup>3</sup>.

## 6. Use

This sensor operates with an interface in standalone mode (without computer), or in connected mode (with computer) after being connected to the interface.

### 6.1 Connection to the interface

The sensor is connected using the connector on the housing:

- Bring the sensor in line with the selected measurement channel
- Insert the sensor into the interface guide system
- Follow the interface guide system.
- Push all the way in

## 6.2 Type of measurement

The sensor measures the concentration of fine dust particles (PM1, PM2.5 and PM 10) in a range from 0 to 500  $\mu\text{g}/\text{m}^3$ .

## 6.3 Operation in standalone mode

### c. Start up

- Connect the sensor
- Switch the interface on

### d. Display

An instantaneous concentration of fine dust particles value is displayed directly on the display screen present on the housing.

## 6.4 In connected mode

### c) Start up

- Connect the sensor to the interface
- Switch the interface on
- Connect the interface to the computer

### d) Display

Measurements are displayed directly on the computer screen, either during an acquisition or using the "display" tab

## 7. Technical characteristics

Channel	3
Particle size	PM1, PM2.5 and PM10
Measurement range	0 to 500 $\mu\text{g}/\text{m}^3$
Resolution	1 $\mu\text{g}/\text{m}^3$
Operating temperature	-10 °C to 60 °C
Stability time	30 s after power-on

## 8. After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

**Free service\***

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux