Ex.A.O.

Capteur Tooxy/FoxyAirNeXT

Data logging

Sensor

Réf :
482 203

Français – p 1

English – p 2

Version : 8110

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1. Généralités

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0 825 563 563*

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1. General

This sensor measures the pH of a solution using a pH probe with a BNC connector (not supplied).


2. Use

This sensor operates with the Tooxy, Foxy and AirNeXT console in standalone mode (without computer) or in connected mode (with computer).

2.1 Connection of the sensor to the console

- Bring the sensor in front of the console connector
- Slide the sensor along the console guide system.
- Push all the way in until it « clicks into place ».

2.2 Start up

Adjustment / calibration :

This action is only possible in standalone mode (USB cable disconnected)

The calibration is made at two points and therefore requires the use of two standard reference solutions, for example pH=4 and pH=7.

Calibration in 2 points

- 1) Do a long press on the selection button: the Solution 1 segment flashes in the violet zone of the display screen.
- 2) Immerse the probe into the first standard solution (for example pH=4)
- 3) Use the potentiometer to adjust the value of the pH (for example to 4)
- 4) Do a short press on the selection button: the solution 2 segment flashes in the violet zone of the display screen
- 5) Rinse the probe and then immerse it into standard solution 2 (for example pH=7)
- 6) Use the potentiometer to adjust the value of the pH (pH=7 in this example)
- 7) Do a short final press on the selection button: the probe OK segment then lights up to show that the probe is working normally.

2.3 Use in standalone mode and in connected mode

In standalone mode, the sensor displays the value of pH with a display precision to 1:10. Indication of the adjustment phase by readouts.

When connected to the computer, the sensor display screen goes off. The measured values are only visible on the computer.

3. Additional equipment

Any pH probe with a BNC connector can be used with this sensor.

4. Technical characteristics

Range	0 to 14 pH
Sensor precision	($\pm 1\% + 2$ LSB) of the value. The final precision depends on the probe.
Measurement time with onboard software	1s

5. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops. For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux