Mécanique

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Mechanics

Acceleration

Réf :
488 060

Français – p 1

English – p 9

Version : 2206

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09 69 32 02 10 *

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1. Description

This wireless accelerometer Jeulin (ref 488 060) is a multi-sensor device detecting acceleration, angular velocity and angle. It provides 6 axis of measurement with a 3-axis accelerometer and a 3-axis gyroscope. An algorithm will calculate the three-axis angle accurately.

The accelerometer is provided with a Bluetooth key and it is compatible with WINDOWS 7, 10 and 11.

Figure 1 : Accelerometer and accessories

Accelerometer et accessories

- 1: Transport case
- 2: Accelerometer wireless
- 3: Bluetooth adapter
- 4: Adhesive strap
- 5: Rechargeable cable
- 6: Adhesive pastilles
- 7: QuickStart Instructions

2. Utilisation

The accelerometer has two main functional modes:

- Wireless mode (via Bluetooth)
- Battery charging mode (via USB-C)

Download software, driver and instructions by using the link below:
www.jeulin.com/installweb, by indicating your activation key provided with QuickStart.

To discover a using video, please scan the QR Code:

2.1 Wireless Mode (via Bluetooth)

If you use the windows 7, you have to install the driver by using the link below with your activation key: www.jeulin.com/installweb

If you use the windows 10/11, the driver is installed automatically.

Insert Bluetooth adapter into the USB port, turn ON the accelerometer.

Figure 2 : Bluetooth adapter connected

Figure 3 : Accelerometer turn ON

If the accelerometer is paired with Bluetooth, the blue LED of sensor will flash and the Bluetooth adapter LED will remain still.

2.2 Pairing

In order to be able to use several accelerometers in the same room, the Bluetooth adapter is paired with the accelerometer. Each of the 2 is provided with a label with its identifier. Be sure to use the associated accelerometer/adapter pairs.

2.3 Battery charging mode

Up to **10 hours** battery life. Rechargeable lithium-ion battery.

You have to connect accelerometer to computer with a USB-C cable in order to charge the battery.

Figure 4 : Battery charging

When the accelerometer is charging, the red LED will remain still, the light will switch off once the charge is completed.

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If the Bluetooth adapter was detected, the software interface would display (Figure 5). Else, a warning screen would display (Figure 6), in this situation, verify that the adapter is functional and retry to connect the adapter.

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
Figure 6 : Warning message

3. Example of experiment

In this part, we will show you an example of experiment which consist to reveal uniformly accelerated rectilinear motion of a car.

3.1 Rectilinear acceleration


Figure 7 : Experiment drawing

In order to study uniformly accelerated motion, we have to stick the accelerometer on the car with adhesive pastille. Then put the car on the test bench, the car is attached by a rope over a pulley to a mass of m which hangs vertically. The plank will be used for stopping the car. Finally, we have to configure the software and launch the acquisition.

3.2 Configuration and acquisition

- Run software "Accéléromètre.exe"
- Drag and drop Y-acceleration **aY** on ordinate axis, and time on abscissa axis. Here, the horizontal acceleration is oriented along the Y-axis (figure 8).


Figure 8 : Software interface

- Define time configuration by clicking on . For this experiment, the time is set up to 5 seconds and 101 points (figure 9).
- Click on **aY**, then "Affichage", multiple calibration exist: $\pm 2g$, $\pm 8g$, $\pm 16g$, select $\pm 2g$ for a better visualization (figure 10).

- Put the car on the test bench, and calibrate accelerometer value to zero:

aY → "Réglage" → "Régler" (figure 11)

Finally, launch the acquisition and release mass m.


Figure 9 : Time configuration


Figure 10 : Choose a suitable display


Figure 11 : Accelerometer calibration

3.3 Results and interpretation

Applying Newton's second law on the car and ignoring friction of car on the bench. We have acceleration:

$$\|a\| = \left(\frac{m}{m + M} \right) * g$$

Right-click on the graph then click on "Droite", place line on top of area which the acceleration is "constant" (figure 11). We have $a = -0,078$, therefore $\|a_{exp}\| = 0,078g$.

Numerical values: $m = 16,1\text{-gram}$ $M = 155\text{-gram}$ $g = 9,8\text{m/s}^2$

	Experimental (Unit: g)	Theoretical (Unit: g)
Acceleration	0,078	0,094

Table : Theoretical and experimental results


Figure 11 : Experimental result


The results are consistent because we neglect the friction. Besides, the accelerometer has an accuracy of $\pm 0,01g$.

Discussing the sources of error may lead to uncertainty of measurement calculations, particularly in connection with the measuring device.

4. Others examples

Different experiments are available in the form of explained and detailed practical work (*pdf). You can find them by scanning the QR codes below.

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The goal is to study linear acceleration on a mechanical bench in order to study the laws of dynamics (force, acceleration, speed, linear movement, etc.).

Voir le TP

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The goal is to study the tangential and radial accelerations of a harmonic oscillator (the heavy pendulum).

Voir le TP

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The goal is to study the acceleration of a mass on an inclined plane in order to study the laws of dynamics, in particular the acceleration of gravity.

Voir le TP

● **TP maquette de rotation**

The goal is to study the moments of inertia as well as the acceleration of a rotating mass around a fixed axis.

Voir le TP

5. Technical Specifications

Battery life	Complete charge	10h
Battery	Lithium ion	See below for more details
Size (LxWxH)	Accelerometer	51,3 mm x 36 mm x 15 mm
	Transport case	140 mm x 100 mm x 37 mm
Measurement Range	Accelerometer	3 axis (x, y, z) $\pm 2g - \pm 8g - \pm 16g$
	Gyroscope	3 axis (x, y, z) $\pm 2000^\circ/s$
	Angle	3 axis (x, y, z) x, z : $\pm 180^\circ$ y : $\pm 90^\circ$
Accuracy/Stability	Accelerometer	Accuracy : 0,01g Stability: 0,005g
	Gyroscope	Stability: 0,05°/s
	Angle	Accuracy x, y: 0,05°

Risks and precautions of Li-ion battery

- DO NOT short circuit (risk of fire and explosion)
- DO NOT overcharge or over-discharge (risk of fire and explosion)
- DO NOT expose to high temperatures
- Batteries are fragile mechanically and electrically
- Only use the charging cord that came with the device
- Avoid charging your device near flammable materials and liquid
- If you notice an overheating, shutoff the charging process
- Avoid charging Li-Ion batteries in series
- Store batteries in a secure, dry and cool place
- DO NOT put batteries on a current lead surface
- DO NOT pierce batteries
- Keep a close watch on batteries during the charge
- A battery could catch fire 10 minutes after a shock or short circuit, keep a close watch on it.
- In case of contact with electrolyte on the skin, wash immediately with soap and water.
- In case of contact with eyes, flush eyes with plenty of cold water and consult an emergency physician.
- Batteries have to be completely discharged before sending them to the recycling

6. After-sales service

The device is under a 2-year guarantee.

For any repairs, contact the **Technical Support** on **+33 (0)2 32 29 40 23**.

The equipment should be returned in our workshops and for any repairs or spares parts, please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

+33 (0)2 32 29 40 23