Ex.A.O

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Data logging

REDY Wireless Sensor

Réf :
488 063

Français – p 1

English – p 14

Version : 4209

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REDY GSR Sensor

Table des matières

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The screenshot shows the 'Acquisition' window with the 'Données' tab selected. Below the graph, there is a section titled 'Temps' (Time) with a 'Synchronisation' sub-tab. It contains several input fields: 'Durée' (Duration) set to 1 s, 'Nb de points' (Number of points) set to 501, and 'Période' (Period) set to 2 ms. There is also a checkbox for 'Acquisition continue' which is unchecked. At the bottom, there are three buttons: 'Lancer' (Launch), 'Arrêter' (Stop), and 'Repère' (Reference). The 'Lancer' button is circled in red.

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Table of Contents

1. Content	13
2. Description.....	13
2.1 REDY sensor	13
2.2 Pulley	14
3. First Use	14
4. Starting Up.....	14
5. Using the Sensor in Connected Mode to a Windows Computer	15
5.1 Installing the software Redy for Windows	15
5.2 Connecting the Sensor to the Software.....	15
5.3 Mechanical calibration of the pulley	Erreur ! Signet non défini.
5.4 Carrying Out an Acquisition Depending on Time	16
6. Using the Sensor in Bluetooth Mode on Tablet	19
6.1 Installing the application Redy for Tablet.....	19
6.2 Connecting the Sensor to the Tablet via Bluetooth	19
6.3 Setting the Acquisition Depending on Time	21
6.3.1 Setting the Intensity Measurement.....	21
6.3.2 Setting the Time Factor	21
7. Technical Specifications.....	23

1. Content

- 1 REDY GSR (Galvanic Skin Response) sensor
- 2 electrodes attached to the sensor
- 1 type C USB cable
- 1 protective Cover

2. Description

2.1 REDY sensor

2.2 Electrodes

3. First Use

When receiving your sensor, and before the first use, it is recommended to fully charge the battery.

Use the USB-C cable supplied to connect the sensor to a USB port of a computer, or a power outlet via a universal charger (like the USB mains charger ref. 805819, not supplied)

Before each practical work, plan to fully charge the device.
The battery life varies depending on the conditions of use and the number of measurements and is designed to ensure a full practical session.

4. Starting Up

Press the Power button.
The device turns on. Place the electrodes on the index and middle fingers of the same hand. The value measured appears on the screen.
Your device is ready for use.

5. Using the Sensor in Connected Mode to a Windows Computer

5.1 Installing the software Redy for Windows


The software Redy for windows is available for downloading on the Jeulin website.

5.2 Connecting the Sensor to the Software

- Connect the sensor to a USB port of your computer.
- Turn on the sensor by pressing on the position I, of the Power button.
- Launch the software Redy for Windows by double-clicking on the icon present on the desktop of the computer.


The sensor is automatically detected. Its symbol appears in the software list of sensors:


5.3 Carrying Out an Acquisition Depending on Time


To carry out an acquisition with the sensor, drag and drop the icon assigned to the measurement (GSR), in ordinate using the mouse:


All the sensor relative indications then appear under the graph:

By clicking on the tab Measure, the value measured by the sensor appears on the screen:

It remains to set the time, by dragging and dropping the timer icon in abscissa in the graph, as shown below:

All the indications relative to setting the time acquisition then appear under the graph:

The screenshot shows the 'Acquisition' window with the configuration panel open. The panel is titled 'Time' and contains the following settings:

- Time function:** Synchronization
- Duration:** 1 s
- Nb of points:** 501
- Period:** 2 ms
- ☐ Continuous acquisition

At the bottom of the panel, there are three buttons: 'Start', 'Stop', and 'Marker'.

Hence, the user can, if he likes, modify the acquisition duration and/or the number of points.

Once the configurations are done, Place the electrodes on the index and middle fingers of the same hand and start the acquisition by pressing on the button « Start » (see above).


6. Using the Sensor in Bluetooth Mode on Tablet

6.1 Installing the application Redy for Tablet


The application Redy is free and available on Google Play or the App Store.

6.2 Connecting the Sensor to the Tablet via Bluetooth

- Turn on the sensor by pressing on the position I, on the Power button.
- On your tablet, launch the application Redy, by pressing on the icon:


- The application is launched, you are asked to connect a sensor:

- Press on yes.

The sensors management space opens automatically in the application.


The sensor appears in the sensors list.

In order to make sure that the desired sensor is detected, you can verify that the address number written on the sensor screen, and that appearing in the application sensors list, match.

This address number is unique for each sensor.


In the example above, it is the "000 - GSR".

- Start the connection by pressing on the desired sensor in the list. A temporary icon appears, then the application confirms that the sensor is connected when the indications below show on the screen:

Note the apparition of the Bluetooth logo next  to the sensor name, as well as the message « Connected » in the end.

- Go back to the main screen by pressing on the back arrow:


6.3 Setting the Acquisition Depending on Time

6.3.1 Setting the Intensity Measurement


The sensor activation, detailed in the previous chapter, selects the sensor automatically for an acquisition. Therefore, no further additional configuration action needs to be taken.

6.3.2 Setting the Time Factor

Click on the icon Setting of the main screen:

Then in the window Acquisition, select Time:

← Acquisition : Time

Duration Synchronization

Duration 10

Unit Second ▾

Continuous acquisition ☐

Allows the acquisition to loop to infinity, this option is useful for short measurements (measurement of a sound for example).

The user can, if he likes, modify the duration and the time unit.

Go back to the main screen by pressing on the back arrow:

← Acquisition : Time

Duration Synchronization


Duration 10

Unit Second ▾

Continuous acquisition ☐

Allows the acquisition to loop to infinity, this option is useful for short measurements (measurement of a sound for example).

Make sure the electrodes are correctly placed and start the acquisition by pressing on the Start button:

7. Technical Specifications

- Skin conductance measurement (GSR Galvanic Skin Response)
- Measuring range: 0 – 10 μS
- Resolution 1 nS
- Backlit graphic screen
- Rechargeable battery Lithium – ion
- USB-C connection
- Bluetooth connection 4.0

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux