Génétique

Maquette

Genetics

Model

Réf :
512 068


Français – p 1

English – p 5

Version : 8007

**Ensemble chromatides
chromosomes**

Chromatides chromosome kit

1 Description

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
1 Description

Single chromosomes formed from a chromatide appear as different colors and sizes of elbow-shaped magnetic sticks. Chromosomes are invisibly magnetized at the centromere and can thus be associated. Their relief and size are adapted to enable easy manipulation. The "elbow" makes it possible to associate two identical chromatides with each other so as to represent double chromosomes (with two chromatides). The magnetic sticks can be identified with dry erasable felt pens and can thus be used to mark genes and alleles differently depending on the activity suggested by the teacher.

2 Composition

14 elements supplied with instructions and support material designed to be easily photocopied:

- 8 "salmon" colored "autosomes" : 4 "long chromosomes" and 4 "short chromosomes" with "elbows" located at distinct locations.
- 4 large green "X sex chromosomes".
- 2 small green "Y sex chromosomes".
- 2 "cellular contours" composed of two flexible rubber tubes provided with magnets at their ends.

Providing many elements makes it possible to work independently, in pairs or in small groups.

3 Complementary equipment

Dry erasable felt pens.
"Demonstration" chromosome models.

4 Technical characteristics

The two faces can be rewritten so that students can work in pairs, each having one side and each having a felt pen color.

5 Usage precautions

This product generates a magnetic field, remember to protect all elements that could be damaged by this magnetic field.

6 Assembly and installation

Models are ready for use and do not need any assembly.

Distribute a set of chromosomes to students or ask students to select the models necessary to solve a specific problem, depending on the proposed activity.

Chromosomes may be placed on views of previously drawn cells (see photocopyable support material provided), and then they can be placed on a sheet and their contour can be plotted so as to keep a written trace of the activity in the file.

The cell wall itself is represented by a flexible tube capable of simulating division of the cytoplasm.

Two sets of 16 models can be combined to form a "fixed" representation of the different allelic combinations, while one model set will be used for a dynamic representation.

7 Use – First implementation


Make photocopies of the support material if you wish (see appendix).

8 Examples of experimental applications, operating methods

8.1 Creating the concept of gene and allele:

Example 1:

4 single chromosomes, two small and two large on which several genes are already annotated, can be presented to students and the students can be asked to associate them in pairs.

Example 2:

The existence and expression of alleles can be demonstrated: for example for the ABO system, manipulation of chromosomes marked differently depending on the alleles provides a means of correlating the existence of several aspects of a character (blood groups A, B, AB and O) with the existence of several alleles (A, B and O) for a gene. Dominance or non-dominance of alleles can also be demonstrated.

8.2 Identification of cellular multiplication steps

Steps of cellular multiplication can be demonstrated starting from a cell containing two pairs of chromosomes (two long and two short). This demonstration may be accompanied by the display of a videogram or an animation so as to make the cellular multiplication steps dynamic.

The 8 chromosomes (4 long and 4 short) may be distributed with a sheet containing a diagram of the initial situation on which the students place 4 single chromosomes at determined locations. Students must then associate them at the centromere so as to form double chromosomes and then separate them into two identical stocks. They can verify that the original number of chromosomes is present in each daughter cell.

The manipulation can be requested "empty" once, with no annotated gene, and then a second time so as to represent a gene by a pair of chromosomes.

8.3 Impact of sexual reproduction on the formation of different and unique human beings

The random distribution of chromosomes carrying different genes and alleles during the formation of reproduction cells can be simulated.

This can be done by monitoring the behavior of two pairs of chromosomes, one of which carries the same two alleles of a gene and the other carries two different alleles.

Example: Rh⁺/Rh⁺ and A/B

Fertilization that causes the formation of genetically different egg cells can also be modeled. The "green" X and Y sex chromosomes, in relative proportions respecting the real size that students can observe in electron microscope photography, can be used to model the determinism of sex as a function of the chromosome carried by the spermatozoon. The support material shown diagrammatically provides a means of placing father and mother chromosomes in the top circles, and then monitoring their distribution in a reproduction cell of each, and finally fertilization by putting the two chromosome contents of the gametes together.

These models can also be used for the following concepts:

- Show that each chromosome contains a DNA molecule with a large number of genes. Each model is represented as being a DNA molecule.
- Illustrate the fact that each daughter cell output from mitosis contains the same genetic heritage as the initial cell.
- Illustrate meiosis and disturbances to it that cause anomalies.

9 Maintenance – Storage

Clean the sticks to remove all traces of marker after use and before storage.
Do not allow felt pen marking to dry, otherwise you might not be able to erase it.

Store protected from the sun and heat.
Do not store behind a window.

10 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries