Génétique

Maquette

Genetics

Model

Réf :
512 070

Français – p 1

English – p 5

Version : 8007

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
1 Description

Chromosomes formed from a chromatide appear as different colors and sizes of elbow-shaped magnetic sticks. They are magnetized and can adhere to or be moved about on a board. The "elbow" is used to put the chromatides into pairs.

The magnetic sticks are coated with a varnish on which erasable dry felt pens can write repeatedly, so that they can be used for permanently marking genes and alleles, depending on the activity presented by the teacher.

Single chromosomes formed from a chromatide are represented by elbowed magnetic supports with different colors and sizes. Chromosomes are magnetized and may be mirrored so as to form double chromosomes. Their relief and size are adapted to enable easy manipulation at the board. The "elbow" makes it possible to associate two identical chromatides with each other so as to represent double chromosomes (with two chromatides). The magnetic sticks are coated with a varnish on which erasable dry felt pens can write repeatedly and can be used for marking genes and alleles, depending on the activity presented by the teacher.

2 Composition

16 "demonstration" models delivered with a manual:

- 8 "salmon" colored "autosomes" : 4 "long chromosomes" and 4 "short chromosomes"
- 4 green "X sex chromosomes"
- 4 green "Y sex chromosomes"

3 Complementary equipment


Dry erasable felt pens.

"Student" chromosome models of the same color and small size for manipulation by students.

4 Technical characteristics

The models are large and very thin : they are designed for easy manipulation and storage, and are easily to see by all students in the class.

5 Usage precaution / safety / regulation

This product generates a magnetic field, remember to protect all elements that could be damaged by this magnetic field.

6 Assembly and installation

Models are ready for use.

Chromosomes can be placed in pairs or singly on the board or they can be placed alone to represent "single" chromosomes, or they can be placed in pairs at the "elbow" that materializes centromeres to represent "double chromosomes", depending on the required activity.

Chromosomes can be placed on views of cells previously drawn on the board and can be moved by the teacher or by a student.

Cell contours can also be displayed on a white board and / or an interactive board using an overhead projector, and chromosomes can be placed on it.

The use of two batches of elements enables "fixed" demonstrations that are more easily understood by a student, with successive stages that remain on the board.


7 Examples of experimental activities

These demonstration chromosomes may be manipulated either by the teacher on the board, or as a complement to the use of student models. In this case, the colors and proportions are the same in the two models, so students can work at their places with their own models, and the activity will be corrected using demonstration models on the board.

7.1 Creating the concept of gene and allele

Example 1

4 single chromosomes, two small and two large on which several genes are already annotated, can be presented to students and the students can be asked to associate them in pairs.

Example 2

The existence and expression of alleles can be demonstrated: for example for the ABO system, manipulation of chromosomes marked differently depending on the alleles provides a means of correlating the existence of several aspects of a character (blood groups A, B, AB and O) with the existence of several alleles (A, B and O) for a gene. Dominance or non-dominance of alleles can also be demonstrated.

7.2 Identification of cellular multiplication steps

Steps of cellular multiplication can be demonstrated starting from a cell containing two pairs of chromosomes (two long and two short). This demonstration may be accompanied by the display of a videogram or an animation so as to make the cellular multiplication steps dynamic.

Four out of the eight chromosomes (2 long and 2 short) can be placed on the board and students must then associate them with the other 4 so as to form double chromosomes and then separate them into two identical stocks. They can verify that the original number of chromosomes is present in each daughter cell.

The manipulation can be requested "empty" once, with no annotated gene, and then a second time so as to represent a gene by a pair of chromosomes.

7.3 Impact of sexual reproduction on the formation of different and unique human beings

The random distribution of chromosomes carrying different genes and alleles during the formation of reproduction cells can be simulated.

This can be done by monitoring the behavior of two pairs of chromosomes, one of which carries the same two alleles of a gene and the other carries two different alleles (example : Rh⁺/Rh⁺ and A/B).

Fertilization that causes the formation of genetically different egg cells can also be modeled. The "green" X and Y sex chromosomes, in relative proportions respecting the real size that students can observe in electron microscope photography, can be used to model the determinism of sex as a function of the chromosome carried by the spermatozoon. The support material shown diagrammatically provides a means of placing father and mother chromosomes in the top circles, and then monitoring their distribution in a reproduction cell of each, and finally fertilization by putting the two chromosome contents of the gametes together.

These models can also be used for the following concepts:

- Show that each chromosome contains a DNA molecule with a large number of genes. Each model is represented as being a DNA molecule.
- Illustrate the fact that each daughter cell output from mitosis contains the same genetic heritage as the initial cell.
- Illustrate meiosis and disturbances to it that cause anomalies.

8 Maintenance and storage

Clean the sticks to remove all traces of marker after use and before storage.
Do not allow felt pen marking to dry, otherwise you might not be able to erase it.

Store protected from the sun and heat.
Do not store behind a window.

9 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries