Génétique

Maquette

Genetic

Scale model

Réf :
512 072

Modèle de mitose
Model of mitosis

Français – p 1

English – p 5

Version : 8005

1 Description

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NOTES

1 Description

Model of mitosis, presented in the form of a cell in bas-relief.

Eight stages of cell division are represented in relief and painted.

Elements represented:

- Genetic material,
- Cellular organelles
- Achromatic spindle
- Centrosomes

End of interphase


Prophase

Prophase

Metaphase

Anaphase

Anaphase

Telophase

2 daughter cells

1.1 The stages of mitosis

During the interphase, the chromosomes are presented in the form of long filaments (chromatin) which cannot be seen individually. The nucleus has a diffuse appearance except for the nucleoles which appear darker. When preparing for mitosis, the genetic material is duplicated, the centrosome (organization centre for the microtubules), located close to the nucleus envelope, also divides.

4 stages can be identified during mitosis: prophase, metaphase, anaphase and telophase.

1.1.1 Prophase

During this phase, the genetic material, (DNA), is condensed into individual structures called chromosomes. The nucleoles disaggregate. Insofar as there are already two identical copies of the genotype in each cell, the chromosomes are constituted of two sister chromatids which each carry the same genetic information. The chromatids of the same pair are held together by a special region of the chromosome called the centromere.

The second major organelle of the prophase is the centrosome. Each centrosome incorporates two centrioles. As for the chromosomes, the centrosome duplicates itself before the start of the prophase. The centrosome which then contains 4 centrioles divides into two and each of the two centrosomes migrates towards one pole of the cell.

1.1.2 Metaphase

The cytoskeleton of microtubules reorganises to form the mitotic spindle, a bipolar structure which stretches between the two centrosomes and is constituted of microtubule fibres.

As soon as the spindle starts to form, each chromosome links up to it by attaching its centromere to several fibres of the spindle. The structure associated with the centromere and forming the bond with the spindle is called the kinetochore.

The assembly of condensed chromosomes can be observed at the equator of the cell to form the equatorial plate. The tensions exercised on each of the kinetochores of a chromosome is gradually balanced out and they align in a plane midway between the two poles. It can be observed the chromosomes are aligned according to their centromere.

1.1.3 Anaphase

The anaphase is a rapid stage of mitosis when the chromatids separate and migrate towards opposite poles of the cell. The centromeres are detached and the chromatids, which then constitute an individual chromosome, each move towards a pole of the cell. During this phase, the sister chromatids separate. They are then "pulled" by the microtubules in the direction of the pole to which they are attached. At the end of the anaphase, the chromosomes are placed in two groups, one at each pole of the spindle.

1.1.4 Telophase

During this period:

- The sister chromatids begin to de-condense and lose their individual appearance.
- The nuclear envelope begins to reform around each group of chromosomes and around the nucleoles.

The spindle gradually disappears. The kinetochore microtubules disappear.

2 Technical characteristics

The cell units are fixed on a metal support

- Format of each support (L x B x H) : 180 x 180 x 30 mm (approximately)
- Total weight: 3 kg (approximately)

3 Use

Manipulating the scale model demonstrates the life of the chromosomes during cell division and makes the students conduct research and reflect.

In an initial stage the scale model can be used to prepare an inventory of the various concepts being demonstrated regarding mitosis.

4 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
0825 563 563

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Service gratuit * :
+ 33 (0)2 32 29 40 50

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Rubrique FAQ

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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