# Immunologie

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## *Immunology*

*Magnetic model*

Réf :  
512 109

Français – p 1

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**0 825 563 563\***

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## 1. Presentation

This model gives an insight into the characteristics of specific humoral immune reaction. Magnetized parts are used to represent the contributors to an immune reaction, and we can then display how B lymphocytes recognize and then neutralize bacteria following contamination through the antibodies that they produce.

The organism continuously recognizes the presence of elements foreign to the organism such as bacteria by means of its immune system.

In most cases following a contamination, phagocytes (category of leucocytes capable of phagocytosis) quickly and efficiently eliminate all pathogenic agents. Otherwise, other slower reactions occur enabling specific recognition of the micro-organism responsible for the beginning of infection.

B and T lymphocytes are involved in specific recognition of antigens (molecules carried by the micro-organism and recognized as being foreign to the organism).

B lymphocytes contain antibodies specific to a single type of antigen. They are capable of recognizing antigens carried or released by micro-organisms.

If such recognition takes place, the B lymphocyte is activated and cloned and will then become capable of releasing large quantities of antibodies with the same specificity. Each antibody has two identical recognition sites. Therefore they are capable of neutralizing the pathogen by forming an immune complex that can subsequently be phagocytosed and eliminated.

A person is seropositive with respect to an antigen when he or she possesses antibodies specific to it.

## 2. Composition

This kit is composed of 30 parts.

## 3. Complementary equipment

The teacher can use the « immunology model » Ref. 512 114 in parallel.

## 4. Technical characteristics

The parts of the model are magnetic. They can thus be placed on a white board and enable interaction with the purpose of displaying some characteristics of the immune defenses of our organism.

## 5. Use - First implementation

The teacher will put the different elements of the model into place on a magnetic white board.

He can add them one by one to support his explanations, or he can put them all into position for a test or an exercise.

## 6. Usage precautions

This product generates a magnetic field, remember to protect all elements that could be damaged by this induced field.

## 7. Examples experimental applications

The purpose of this model is to materialize the effects of B lymphocytes, the antigen antibody specificity and to demonstrate the formation of antigen-antibody complexes resulting from them.

Manipulations may be done by the teacher alone or by students called to the blackboard. For example, they can build up a summary diagram at the end of a chapter.

### Specificity of antigens

The model presents two types of bacteria that have a different antigen on their surface. Thus, a different antigen will be found in the organism depending on the microorganism responsible for the infection.

Required equipment: the two types of bacteria.

### Activation of effective B lymphocytes

There are two types of B lymphocytes in the model. Only one of the two will be used in reactions against the infection. The only one that will be activated will be the lymphocyte that has the specific membrane antibody of the antigen carried by the pathogen responsible for the infection. This is done by putting the two types of lymphocytes and one bacteria side by side. The purpose of this operation is to demonstrate complementarity between the antigen carried by the bacteria and the membrane antibody carried by one of the two lymphocytes.

Required equipment: two lymphocytes with different specificity and one of the two types of bacteria.

### Multiplication of lymphocytes specific to a recognized antigen.

In the model, there are 6 copies of the first type of B lymphocyte and 2 of the other. They can demonstrate that after specific recognition of the foreign antigen, the population of specific B lymphocytes is strongly amplified. It can be demonstrated that starting from an activated B lymphocyte, several clones of this same B lymphocyte are obtained.

Note that this multiplication takes place within some organs and particularly in the lymph nodes.

Required equipment: all B lymphocytes specific to the antibody carried by the bacteria present.

### Production of antibodies specific to the antigen

In the model, there are antibodies associated with each B lymphocyte (9 copies of one and 3 of the other). They demonstrate that once the specific B lymphocytes have been activated, they become capable of producing a large quantity of antibodies with the same specificity, in other words that are capable of recognizing the same antigen, in the blood. This is done by bringing B lymphocytes specific to the antigen together with the antibodies that they are capable of producing in large number.

Required equipment: B lymphocytes and antibodies specific to the antigen.

### Action of antibodies

Each antibody has two identical recognition site. Therefore they are capable of neutralizing the pathogen by forming an immune complex that can be subsequently phagocytosed and eliminated. With the model, bacteria responsible for the infection can be brought together with antibodies specific to them. One or several immune complexes will be formed, due to the complementarity of forms between bacterial antigens and antibodies.

Once the bacteria have been immobilized by the antibodies, they lose their pathogenic capacity and will be quickly eliminated by phagocytosis.

Required equipment: all antibodies specific to bacteria.

#### Seropositivity

When a person has the specific antibodies against a given pathogen, this person is said to be seropositive. The blood content and the presence of specific antibodies can be represented to illustrate this concept.

Required equipment: the two types of antibodies.

## 8. Maintenance - Storage

Magnetized parts must be kept in a plastic envelope protected from light to prevent colors from fading.

## 9. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

**Free service\***

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux