Génétique

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Genetic

Blue LED transilluminator

Réf : 527 004

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Blue LED transilluminator

Français – p 1

English – p 4

Version : 0210

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0 825 563 563*

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1. Description of Blue LED transilluminator

The transilluminator is composed of :

- blue light source 470 nm
- 12V power supply
- orange filter

Dimensions of the transilluminator :
115 x 80 x 160 cm

Dimensions of orange filter :
22 x 9.5 cm


2. Detection of DNA by agarose gel electrophoresis

2.1 Principle

This transilluminator coupled with a fluorescent dye of the Gelgreen® type makes it possible to reveal the presence of DNA bands in agarose gels with a high level of sensitivity.

The device advantageously replaces the conventional visualization systems by azure blue A staining. Easier and faster because it eliminates discoloration baths and more precise thanks to greater sensitivity to DNA.


Gelgreen is a dye, non-toxic to health, specific for DNA. Composed of a large molecule, Gelgreen® does not enter cells, so it is very different from a DNA intercalator like EtBr. Gelgreen® is integrated either before or after migration, it has a high affinity for DNA, so the DNA fragments present in the gel will load.

Under the effect of the blue light from the transilluminator, the dye fluoresces by re-emitting a yellow light.

The role of the orange filter is to "stop" blue light (470 nm) and allow yellow light (approximately 540 nm) to pass.

2.2 Implementation

The implementation just takes a few seconds. After migration, the transilluminator is placed under the gel either directly under the tank, or using its support, the gel is deposited directly on the transilluminator.

We turn on the transilluminator in the dark, or in the darkroom (ref. 527010), the DNA bands appear clearly within a few minutes (3 to 5 minutes).

No UV lamp is needed. The manipulation can be carried out without difficulty in class by the students.

Transilluminator with the DNA electrophoresis tank, reference 591031

3. Transgenesis : visualization of GFP proteins


The transgenesis practice consists of the introduction of a new fragment of genetic material (DNA) into an organism (yeast or bacteria). If this DNA is correctly expressed, then organism will be provided with a new identifiable trait.

This transgenesis uses a small circular DNA, a plasmid, which is naturally present in some yeasts and bacteria. The plasmid was artificially prepared in laboratory and carries GFP gene from jellyfish which codes for a fluorescent protein (GFP = "Green Fluorescent Protein"). So when transgenesis is successful, GFP gene is expressed and microorganism that produces the GFP protein becomes fluorescent when illuminated with a specific wavelength. Biotechnology laboratories have designed numerous GFP plasmid variants for specific uses (CopGFP, daGFP, PEGFP etc.).


The * daGFP protein

This plasmid codes for a fluorescent protein called daGFP, this protein derived from native GFP has its own spectral characteristics. It can be conventionally used in UV excitation (max 395 nm) and also has a second excitation peak at 510 for an emission at 521 nm to 540 nm (green).

Using the 470 nm blue light transilluminator, viewing is facilitated by coupling with an orange filter, which cuts emissions below 510 nm.

Yeast in YPG box → Yeast processed with GFP

Transilluminator and orange filter or darkroom for transilluminator

4. After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact :

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux