Géologie

Geology

**Réf :
545 017**

Français – p 1

English – p 7

Version : 6103

Radiomètre

Radiation meter

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0 825 563 563*

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Earth observation satellites (Spot, Landsat, Météosat) do not usually take photographs of the Earth, but they measure the luminance (L) of the different overflown points using on-board radiation meters. Luminance is measured in several spectral bands for each point. These various measurements are sufficient to differentiate soils and plants on the Earth's surface.

The principle of these measurements can be materialised by moving the radiation meter at a constant altitude over a collection of different objects and plants arranged in line, and observing variations of luminance firstly in one spectral band and then in others during several passes over the same zones.

These measurements obtained point by point as for satellites, model the principle of data acquisition and can demonstrate that:

- it might be possible to differentiate objects that have the same response in one spectral band, in other spectral bands;
- identical objects have identical spectral responses.

The use of the radiation meter is justified particularly in the Final high school year Scientific class, in theme 2 of the Natural & Life Sciences specialty:

« Contemporary global challenges - Atmosphere, hydrosphere, climates: from the past to the future ».

More precisely, the radiation meter is a means of characterising the nature or composition of soils.

1. Principle and description

1.1 Principle

The theoretical study shows that the response of the radiation meter is proportional to luminance of the target along the line of sight, in other words the intensity radiated along this direction per apparent unit area of the target. Luminance is expressed in Watts per steradian and per square meter.

The luminance of a surface varies with the wavelength λ of the radiation. Therefore, we need to define the spectral luminance $L(\lambda)$ for a given wavelength.

It is important to know that, if the surface aimed at by the radiation meter is uniform, the result of the measurement (allowing for precision) is independent of:

- the distance from the radiation meter to the surface,
- the dimensions of the observation field,

Therefore, luminance is a fundamental magnitude in remote detection.


The reflectance R of a surface is the quotient of the reflected radiating power and the received radiating power. Therefore, it is between 0 and 1 but it is often expressed as a percentage. Unlike luminance, reflectance is a magnitude that is characteristic of the environment. It is the basis for recognition of substances. Radiation metering can be used to make a quantitative determination. Spectral reflectance $R(\lambda)$ is the reflectance defined for each wavelength.

With a constant illumination of the target, the following are compared for a given spectral band and for the same line of sight:

- the luminance L of the substance studied,
- the luminance L_0 of a matt white surface (= white reference) with a known reflectance (for example, 180 g/m² of matt white Canson paper with reflectance of about 0.92 for the wavelength range from 400 to 1100 nm).

The JEULIN radiation meter is capable of measuring the reflectance of any object towards which it is aimed.

2. First start up

Radiation meter

- Start up: 12V power supply module provided (reference : 281 243).
- The **Min** and **Max** buttons control amplification, and are used to set the measurement range for a reading between 0 and 100% regardless of the spectral band chosen.
 - If the signal is too strong when the « min » or « max » is being set, then the letters Hi appear on the screen showing that it is impossible to make the measurement.
 - If the signal is too weak when the « min » or « max » is being set, then the letters Lo appear on the screen showing that it is impossible to make the measurement.
- The two half disks turn independently so that the user can choose the working spectral range (R for Red, G for Green, B for Blue, IR for Infra Red and NF for No Filter) and the aperture (more or less large).

3. Parameter settings before measurement

3.1 Calibrating the radiation meter

3.1.1 Setting the minimum: 0%

- Select the spectral band by turning the filter set (NF, IR, R, G, B)
- Put a piece of matt black 160g/m² Canson paper under the radiation meter that is switched on (refer to the zero setting image)
- Put your light source into place and switch it on.
- Keep the **Min** button pressed for a few seconds to obtain the 0% display (within 10%).

Setting the zero

3.1.2 Setting the maximum: 100%

Keep the previously chosen spectral band.
Hold the matt black Canson paper in place.
The source must not move from the previous position.

Place the substance to be studied on 160 g/m² matt black Canson paper and under the radiation meter, the substance must be placed in a plane parallel to the radiation meter.

Cover the target with a sheet of 180 g/m² matt white Canson paper
Keep the « max » button pressed for about 2 seconds until the value of about 100% appears.

Setting the maximum

3 cases can arise:

- Case No.1:

The displayed value is close to 100% (it can vary by plus or minus 8% depending on the stability of the lighting), the setting is terminated for the studied spectral band.

- Case No. 2:

The display shows **Lo** which means that the light power is badly adapted. The diaphragm has to be changed. In this case, the zero setting must be made again and then the maximum setting must be made again.

- **Case No.3:**

The display shows **Hi** which means that the light power is badly adapted. The diaphragm has to be changed. In this case, the zero setting must be made again and then the maximum setting must be made again.

3.2 Measurement

Remove the white sheet and read the reflectance value on the screen.

This value as a percentage is the reflectance of the substance studied in the spectral band corresponding to the filter used.

Repeat these operations with the other available filters to obtain 4 points on the reflectance curve of the substance. An example is given below.

3.3 External analogue outputs

The analogue output signal is 0-1 V for a connection to a CAT console or to a multimeter.


To obtain the reflectance, record the voltage in Volts and convert it to %.

For example, if the voltage is 0.2 V, the reflectance is 20%.

It is essential to calibrate the instrument for each wavelength.

The wavelengths for each filter are given at the end of these instructions.

Example acquisition of the reflectance of earth placed in a Petri dish as a function of the wavelength.

4. Precautions for use

Radiation meter measurements require very strict working conditions:

- The target aimed at must be on the optical axis of the radiation meter and within the observation field.
- Avoid shadows created on the target by the radiation meter, the operator, spectators or surrounding objects.
Changes in the position of the experimenter, his hands or spectators can change illumination conditions.
- Measurements are usually made in varying width spectral bands, rather than in monochromatic light. The measurement result then depends on the spectral distribution of the light source and is slightly different depending on whether the measurement is made in daylight (different at midday and at sunset) or under an incandescent lamp or a fluorescent lamp or even in the light of a room with coloured curtains.
- Since the brightness of the target can vary, specular reflection (as from a mirror) can mask diffuse reflection or transmission if the angle of inclination relative to the light is badly chosen.
- Put matt black Canson paper behind the targets to prevent hindering reflections.

Comments:

- Illumination must be the same on the white reference, the black reference and on the substance studied; therefore when working in daylight, the measurement must be made immediately after removing the white reference (the illumination can vary very quickly when the sky is cloudy).
- If illumination is insufficient, use an incandescent lamp placed close to the radiation meter illuminating the target at an angle of incidence of about 45 degrees.
- It is important to know that reflectance from plants varies depending on whether the measurement is made on a single leaf or on two superposed leaves, on a more or less light, healthy or faded leaf, on a whole plant or on a set of identical plants. The difference is due to the variation of the orientation of the leaves, the presence of shade or possibly soil between the leaves or the plants.

5. Technical characteristics

5.1 Instrument

Dimensions: 110 x 120 x 50 mm
Mass: 230 g

External analogue outputs 

The 0-1 V output displays the measurement result on a voltmeter or on any CAT console.

5.2 Filters

When you want to change the filter, rotate the filter-holder disk to obtain one of the 5 positions. The mark located on the top face of the radiation meter shows which filter is placed under the detector.

The following table gives the correspondence:

MARK	FILTER	(nm)
NF	No filter	Ø
IR	Infrared	800-1100
R	Red	600-700
G	Green	500-600
B	Blue	400-500

With a 5-position diaphragm, a sensitivity adjustment complementary to the radiation meter adjustment can be made.

6. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux