**Méthode ESAO®**

***ESAO® Method***


**Réf :  
554 004**

**Français – p 1**

**English – p 5**

**Version : 7110**

**Stimulateur  
*Stimulator***

## 1. Généralités

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**0 825 563 563\***

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## 1. General information

Manufactured like an Initio case, this stimulator helps give out single, double, multiple or ascending stimulation pulses.

It enables complete adjustment of the stimulation parameters (response time, intensity, duration, interval, number of pulses, variable intensity) while synchronising the acquisition.

## 2. Use

### 2.1 Connection

Power the stimulator using a 12V/500mA power supply (Part no. 281 512), not supplied.

The stimulator has 2 sockets on the left that are used to connect the stimulator to the synchro measurement device.

#### With ESAO:

Connect the red socket of the stimulator to the red socket of the synchro located on top of the console.

Connect the black socket of the stimulator to the black socket (earth) located on top of the console.

On the right side, the stimulator has **two sockets, one yellow and one black that are the outputs**. Connect these outputs to the set up through a nerve shell, nerve desk or stimulating electrodes.

## 2.2 Parameterisation

### ❖ Mode selection

Select the stimulation mode first: single, double, bursts or multiple.  
This selection may be made using the grey button located below. An LED indicates the current selection. The default mode is single pulse.

Every mode provides access to different adjustment parameters in relation to the selected stimulation mode.

#### Single pulse mode:

Adjustment of the response time d, duration of stimulation t, intensity of the stimulation U1.

#### Double pulse mode:

Adjustment of the response time d, duration of stimulation t, interval  $\Delta$ , intensity of the stimulation U1. The interval corresponds to a complete period (stimulation period + non stimulation period).

#### Burst mode:

Adjustment of the response time d, duration of stimulation t, interval  $\Delta$ , intensity of the stimulation U1.

In this mode, the number of constant intensity pulses is determined by the duration of depression of the start button.

#### Multiple mode:

Adjustment of the response time d, duration of the stimulation t, interval  $\Delta$ , number of pulses n, intensity of the first stimulation U1, intensity of the last stimulation U2.

This mode may be used with an increasing, decreasing or fixed intensity of stimulation. In this case the stimulation in bursts can be given with a known number of pulses.

After the mode is selected, it is advisable to adjust the various parameters. Only the parameters available in the selected mode are accessible.

### ❖ Adjustment of parameters

#### Choice of the parameter

Under the screen is a set of 7 LEDs corresponding to the parameters. Use the 2 grey buttons located on each side of the set of LEDs to change from one parameter to the other. The values displayed on the screen correspond to the parameter selected (green LED on).

#### Modifying the value of the parameter

Using the 2 grey buttons (+/-) located on the right of the screen: increase or reduce the value.

#### The adjustment ranges for these parameters are:

Response time	From 0.5 ms to 999 ms
Time	From 0.1 to 999 ms
Interval*	From 0.2 ms to 999 ms and from 1 s to 100 s
Intensity	From 0.1 V to 5 V
Number of pulses	From 1 to 100

\* The duration of the interval must be greater than the time t of the stimulation.

❖ **Special case of the adjustment of the interval** △

If one wishes to set the interval in seconds, it is necessary to go through the adjustment of parameters in ms and increase the value by pressing on the (+) button until the maximum value is obtained. Then the transition to seconds is done automatically.

❖ **Initialising the stimulation**

Press the start button located on the bottom right to launch the stimulation.

In burst mode and only in this mode, the duration of depressing this button corresponds to the duration of the stimulation.

### 3. Technical specifications

- Initio Type Case
- 2-stage adjustment using push-buttons: Choice of the parameter followed by adjustment
- Indication of the parameter selected using LEDs
- Can be used with all systems having a Synchro output
- Display on a 12 mm 3.5-digit LCD screen
- 1 synchro output, 1 signal output
- Power supply through a 12 V/500 mA mains unit **not supplied**
- Dimensions: 10 x 10 x 4.5 cm
- Weight: about 200 g

### 4. Voltage adjustments

- In case of a drift and only in this case, it is possible to readjust the operating voltages.
- Connect the stimulator to a multimeter using the signal output sockets located on the right-hand side.
- While pressing the 2 buttons on the left for the mode selection and parameter selection (arrow towards the left) connect the stimulator to the mains.
- Then, using the parameter change button select U1 (50 mV). Adjust the value using the +/- buttons and validate by pressing start.
- Then, using the parameter change button select U2 (5V). Adjust the value using the +/- buttons and validate by pressing start.

The adjustment is completed.

### 5. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.  
For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
468 rue Jacques Monod  
CS 21900  
27019 EVREUX CEDEX FRANCE

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

### Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux