Physiologie

Physiology

**Ref :
554 009**

Français – p 1

English – p 3

Version : 6010

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2 Notes

- The 6 small wells are to hold limewater, and their perforated lids allow gas exchange to take place.
- The 6 closed large wells are to hold the small wells and the samples to be studied.

The carbon dioxide concentration in the air increases due to the respiration of the study organisms. The heavier carbon dioxide descends through the perforated lid and reacts with the limewater.

The diameter of the wells is adapted to prevent the samples (notably mobile samples such as mealworms, etc.) from climbing down the walls and getting caught between the observer's eye and the limewater.

Examples of organisms that can be placed in the well are: carrots cut into pieces, mushrooms cut into pieces, mealworms, maggots, etc.

3 Instructions

With the pin, perforate each of the lids to the small wells; make around 10 to 12 holes. Fill each small well about half full with limewater and close it with the perforated lid.

- Place each of the small wells prepared in this manner in a large well, and place the sample to be tested on top of the perforated lid. Close the large well with its lid.
- Depending on the sample being studied, allow the respiration to take place for several minutes or several tens of minutes.

4 Storage and care

- Wash the apparatus with a little bit of detergent.
- Clean the limewater wells with slightly acid water (diluted citric acid solution item no. 107 221 or diluted ethanoic acid solution item no. 106 106) to remove the lime residue.
- Do not use a scouring pad, as it will scratch the wells.
- To avoid any traces of lime, rinse with demineralised water, allow to dry and store in a dry place.

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*