Méthode ESAO[®]

ESAO[®] method

**Ref :
554 041**

Français – p 1

English – p 3

Version : 6006

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4. Instructions

Step 1: Sit down comfortably near a table and relax in this position for a few minutes. Raise your sleeve up to the elbow and make sure that your clothes don't restrict your arm.

Step 2: Position the cuff and make sure that it is not too tight or too loose. Place the stethoscope probe on your brachial artery that is slightly on the inside surface of the arm above the crook of the elbow.

Step 3: Place the earpieces of the stethoscope to your ears and close the screw valve by rotating it clockwise. Keep the arm extended with the elbow slightly bent, palm facing up and relaxed.

Step 4: Compress the pump bulb using your free hand until the manometer indicates a value greater than about 30 mm Hg of your expected blood pressure. If you don't know the value, the minimum pressure is usually 200 mm Hg. Stop pumping and rotate the screw valve slowly counter clockwise to reduce the pressure at a rate of 2 to 3 mm Hg per second.

Step 5: As the pressure drops, you will hear a low but clear drumming sound, which progressively increases as the pressure continues to drop, with varying intensity and type until it disappears completely. When the sound has disappeared, open the screw valve to deflate the inflating system rapidly and completely. Watch attentively the manometer of the measuring device while listening to the sound with the stethoscope. Note down the value indicated by the manometer at the point where you begin hearing the drumming sound, as this represents your maximum systolic blood pressure. The diastolic blood pressure (minimum pressure) is determined by the complete disappearance of the sound. Note down the value indicated by the manometer at this instant.

Step 6: Note down these two values as well as other pertinent information, such as the date and time of measurement.

5. Some important points

- ❶ Do not tighten the cuff too much as this makes it difficult to determine the diastolic blood pressure.
- ❷ Do not operate the pump bulb with the hand of the arm wearing the cuff.
- ❸ Do not reduce the pressure too fast. Rotate the screw valve gradually and carefully while releasing air, otherwise you will have an inaccurate pressure value.
- ❹ If you are unable to read the manometer correctly or if you missed noticing the sound at any time during measurement, do not try to reinflate the cuff. In such a case, deflate the cuff completely and wait for 2 or 3 minutes before restarting.

Usually, systolic pressures are between 100 and 120 mm and diastolic pressures between 60 and 80 mm in young and healthy patients. The blood pressure varies according to factors such as age, sex (women usually have a pressure that is lower than men), behaviour, weight, muscular development, anxiety and fatigue, all these having a significant effect on your pressure.

If you have any questions about how to measure your blood pressure, please contact your doctor, who will give you appropriate advice.

NEVER FORGET THAT IT IS ONLY YOUR DOCTOR WHO IS QUALIFIED TO EVALUATE THE SIGNIFICANCE OF YOUR BLOOD PRESSURE MEASUREMENTS.

6. After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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Service gratuit * :
+ 33 (0)2 32 29 40 50

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Rubrique FAQ

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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