Microscopie

Microscopy

**Ref :
571 220**


Français – p 1

English – p 5

Version : 9009

Microscope monoculaire à LED 40x1000

Monocular microscope

Légende :

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B Objectifs

C Platine carrée

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F Tête monoculaire

G Tourelle 4 places

H Surplatine

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J Mise au point fine

K Bouton de réglage de la platine

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
L Light intensity adjustment

1 Technical characteristics

Head	monocular, tilted at 45°, pivoting at 360°, non detachable
Eyepiece	x10 set wide field Inverted 4 position turret
Achromatic objectives	x4 N.A. 0.1 ; x10 N.A. 0.25 ; x40 N.A. 0.65 ; x100 N.A. 0.85
Magnification	x40 to x1000
Stage	Fitted with a 125x125mm slide holder Rack and pinion adjustment with separate micrometric and macrometric commands
Focusing	
Rack stop	Yes
Diaphragm	Iris with filter with filter holder
Blue filter	Yes
Cover	Provided
Condenser	Adjustable height 1.25 Abbe condenser
Lighting	LED lighting, wheel adjustment separate from switch Rechargeable batteries
Power	and separate 5.5V-300mA transformer
Weight	5 kg

2 Assembly

- Remove the microscope from its box and packaging, carefully place it on a stable work bench.
- Remove the plastic bags and the dust covers from the adaptors.
- Place the monocular head in the tube adaptor; tighten the knurled screw with fingers.
- Familiarise yourself with the mechanics of your microscope. Gently handle each element to see how it behaves and what effect it has. Connect the power cord correctly using the power supply provided with the microscope.

Note :

- 1) The microscope must have an earth connection.
- 2) Make sure that the supply voltage is adequate for the nominal voltage of the microscope.

The batteries have an initial charge time of 5 to 6 hours. The standard charge time is of 4 hours for an autonomy of 4 hours.

3 Use and adjustment

- To move the microscope, hold it by the inclined portion of the stand only, placing the other hand underneath.
- Turn on using the switch at the back of the base, turn the light dimmer switch on the right of the base to obtain the desired illumination.

- Place the sample (slide) to be examined delicately on the stage, cover slip towards the objective lens. Delicately pinch the sample (slide) using the slide holder.
- The 360° pivoting observation head and the 45° inclined tube provide an excellent ergonomics and allow prolonged use.
- A light condenser is found under the stage. The intensity of the incident light may be modified when the diaphragm aperture is adjusted.

Note: The diaphragm aperture is not designed to adjust the luminosity; the luminosity is adjusted using the luminosity adjustment knob.

Turn the objective holder to change the 4x or 10x objective, and make sure (by hearing a “click”) the objective is in the light path.

To prevent the objective from touching the sample when adjusting the focus, turn the coarse focus knob until the sample is at approximately 3 mm from the objective. Slowly turn the coarse focus adjustment knob until a clear image is obtained; then use the fine focus knob to make the sample image as clear as possible.

4 Maintenance

- Dust protection: cover the apparatus with its cover after use. If dust deposits itself on the lenses, blow using a syringe or a dust remover (Ref. 574 198) or use a sable hair brush: the lenses are fragile and may be scratched by dust.
- Fingerprints: do not touch the lenses. However, if they have fingerprints, remove them with a soft cloth or an optical cleaning paper (Ref. 803 092).
- Humidity: humidity, especially combined with temperature, is particularly harmful for mounted optics. Microscopes should be stored in a dry place, and in particularly exposed climates, in their polystyrene boxes with a dessicator bag.
- Dismounting: a microscope is a precision apparatus with delicate machining. If the dismounting is necessary (maintenance), return the instrument to the company JEULIN. The lower plate that closes it should only be opened by an electrician as the case must remain protected at the electrical safety level. Only dismount when the microscope is disconnected.

The shelf life of the LED is of 10,000 hours, in case of non-function of the lighting and after verifying that the batteries are charged, send the microscope back to Jeulin.

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)825 563 563

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries