Observation

Microscopes

Observation

Microscopes

Réf :
571 235

Français – p 1

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
+33 (0)825 563 563*

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Caption:

- | | |
|--|---|
| A Wide field eyepiece with index | G Inverted quadruple nosepiece |
| B Achromatic objectives | H Double layers mechanical stage with coaxial adjustment knobs |
| C Square stage | I Coarse focus |
| D Adjustable Abbe condenser | J Fine focus |
| E Illumination | K Dimmer |
| F Monocular head with built-in camera | L Base with batteries |

1 Technical specifications

Head	Monocular, 30° inclined, 360° rotatable, non removable, with built-in 1,3Mp digital camera
Eyepiece	Wide field x10
Achromatic objectives	Inverted quadruple nosepiece x4 N.A. 0.1 ; x10 N.A. 0.25 ; x40 N.A. 0.65 ; x100 N.A. 1.25
Magnification	x40 to x1000
Stage	Double layers mechanical stage 120x120mm
Focusing	Rack and pinion adjustment with separate micrometric and macrometric controls
Rack stop	Yes
Diaphragm	Iris with filter holder
Blue filter	Yes
Cover	Yes
Condenser	Abbe N.A. 1.25
Illumination	LED with separate adjustable brightness
Power supply	Rechargeable batteries and separate 5,5V-300mA transformer
Weight	5 kgs

2 Set-up

- Remove the microscope from its box and packaging, carefully place it on a stable work bench.
- Remove the plastic bags and dust covers.
- Screw the x4, x10, x40 and x100 eyepieces.
- Become familiar with the mechanics of your microscope. Gently handle each element to see how it behaves and what effect it has.

Normally, the microscope is supplied with loaded batteries. However, it is recommended to load them so that the illumination remains optimal.

3 Use and adjustment

- When using the microscope, pick it up only with the inclined arm of the body, and place your other hand under the base.
- Use the switch behind the base of the microscope to turn it on, turn the dimmer located on the base to obtain the desired light intensity.
- Place the prepared slide on the stage, cover slip towards the objective lens.
- The 360° rotatable head and the 30° inclined tube are ergonomically designed for prolonged use.
- A light condenser is found under the stage. The brightness may be modified when the diaphragm aperture is adjusted.

Note: The diaphragm aperture is not designed to adjust the brightness; the brightness is adjusted using the brightness adjustment knob.

- Turn the objective holder to change the 4x or 10x objective, and make sure (by hearing a “click”) the objective is in the light way.
- To prevent the objective touch the prepared slide when focusing, turn the coarse focus knob until the prepared slide should be closed to 3 mm from the objective. Turn slowly the coarse focus adjustment knob until a clear image is obtained:
then use the fine focus knob to make the prepared microscope image as clear as possible.

4 Camera set-up

Insert the CD provided with your videomicroscope into the relevant drive.

Connect the camera to a USB port of your PC.

Follow the USB Windows Standard set-up procedure. At the end, your microscope is recognized.

5 Operating software

Your videomicroscope is supplied with an operating software available on the provided CD.

Your microscope is compatible with the Picturio software.

6 Maintenance

- Dust protection: After use, place the cover over the apparatus. If there are dust deposits on the lenses, blow them away with a syringe or a dust remover (Ref. 574 198) or use a martin hair brush ; the lenses are fragile and may be scratched by dust.
- Fingerprints: Do not touch the lenses. If, however, if they do have fingerprints on them, remove them with a soft cloth or an optical cleaning paper (Ref. 803 092).
- Humidity: Humidity, especially combined with temperature, is particularly harmful for mounted optics. Microscopes should be stored in a dry place, and in particularly exposed climates, in their polystyrene boxes with a dessicator bag.
- Dismounting: A microscope is a finely-tooled apparatus. If it is necessary to take it apart (for maintenance), return the instrument to the JEULIN company or a qualified repair shop. The base plate that seals the microscope may only be opened by an electrician because the body must remain protected according to electrical safety standards. Do not take the microscope apart unless it is unplugged.

The shelf life of the LED is of 10,000 hours, in case of non-function of the lighting and after verifying that the batteries are charged, send the microscope back to Jeulin.

7 After-Sale Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries