Observation

Stéréomicroscopes

Observation

Stereomicroscopes

**Réf :
571 330
571 331**

Français – p 1

English – p 5

Version : 3107

**Loupes zoom autonomes
Zoom stereomicroscopes**

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0825 563 563 *

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
1. Presentation

These high performance zoom stereomicroscopes are provided with a zoom system graduated from 1 to 4. The total magnification varies from x10 to x40. These are meant to observe any biological or geological material, which cannot be observed with a microscope due to its thickness and/or its opacity.

It exists in 2 versions: binocular and trinocular model. They are both stand-alone device with its rechargeable batteries (supplied with). The trinocular model can be used with a digital eyepiece camera so that the whole class can see the display.

2. Description

Trinocular model, ref. 571 331

Caption:

- | | |
|---|---------------------------------------|
| A Eyepiece tubes with dioptic adjustment | F Zoom adjustment |
| B Trinocular head | G Macrometric focusing |
| C Zoom objective | H Illumination switches |
| D LED incident and/or transmitted illumination | I Brightness knob |
| E Video tube | J Base with built-in batteries |

3. Specifications

	Ref. 571 330	Ref. 571 331
Head	Binocular	Trinocular
	45° inclined, 360° rotatable, with eyeshades	
Eyepieces	WF10x, Ø 30,5 mm	
Interpupillar distance	55 - 75 mm	
Dioptric adjustment	Both tubes	
Objectives	Graduated zoom from 1 to 4 Graduations: 1 - 1,5 - 2 - 3 - 4 with a separate knob	
Magnification	x10 – x40	
Field of view	23 mm	
Working distance	61 mm	
Stage	1 white/black and 1 transparent	
Illumination	Both adjustable incident and transmitted LED	
Power supply	Rechargeable batteries (supplied with) and 4,5 V / 600 mA transformer	
Dust cover	Yes	

4. Unpacking and assembling

4.1. Unpacking the stereomicroscope

Carefully open the styrofoam packaging and check that all components have been taken out the packaging. Do not throw out the packaging before identifying all the components.

4.2. Assembling the microscope

When manipulating the components, and particularly optical parts, take care to not touch the lenses with your bare fingers or hands. Any fingerprints or grease spots will reduce the image quality.

After unpacking the support, install it on a stable working surface and remove the various protections items (do not touch the surface of the lens). Fix the head and install the eyeshades on the eyepieces.

With the trinocular model, mount and fix the head. Place the eyeshades onto the eyepieces. Unscrew and remove the objective cover. When pushing the lever on the left of the head, you will work with a binocular model.

It exists 2 ways to supply the stereomicroscope:

- stand-alone with the supplied rechargeable batteries,
- with the transformer.

First use:

Before connecting the instrument and switching it on, make sure that the operating voltage of the stereomicroscope is the same as the power supply voltage.

Batteries specifications:

Loading time = 4 hours ; Stand-alone time = 4 hours.

4.3. Alignment and working

4.3.1. Interpupillary distance

Move the two eyepiece tubes until only one circular field is visible through the two eyepieces. If two circles appear, the interpupillary distance is too large, and if 2 overlapping circles appear, the interpupillary distance is too small.

4.3.2. Focusing the stereomicroscope

Chose the lowest magnification and adjust the image with the macrometric knob until having the clearest image. Increase the magnification step by step and readjust each time the image until having the required image.

4.3.3. Stage removing

In regards of what you will observe, you will prefer a clear or dark stage. With this stereomicroscope is provided 2 stages: the one in transparent frosted glass, the other is a reversible black and white stage. In order to remove this stage, unscrew without removing it the screw in front of the stage. Then, lift the stage in the back.

The transmitted illumination is only used with the frosted glass stage.

4.3.4. Illumination

For your comfort, this stereomicroscope is supplied with separate incident and transmitted illumination. Adjust the brightness with the illumination knob.

4.3.5. Digital eyepiece camera using

With the binocular model: unscrew without removing it, the screw holding the eyepiece, then, insert the digital eyepiece camera tube in place. If an adapter is required, then insert it depends on the digital eyepiece camera instructions. Adapters are supplied with the digital eyepiece camera.

With the trinocular model: screw the metallic tube onto the video tube. Insert the digital eyepiece camera tube. Depending on the digital eyepiece cameras, you can use the adapter supplied with. Take care that the lever on the left of the head is not push, otherwise you will work with a binocular model.

5. Cleaning and maintenance

Protection against dust: cover the instrument with the protective cover provided if it remains unused for a long period. Never leave the eyepiece tube open and unprotected. Either leave the eyepiece in the tube (which is recommended if the microscope is used very frequently) or cover it with wrapping paper or a protective cap. It is recommended that eyepieces and other optical accessories should be kept in a dry box to protect them from dust and humidity when they are not in use.

Protection against humidity: the instrument must be kept at a distance from every water source, pipe and sink. Humidity in the room in which the instrument is to be used must be as low as possible (the relative humidity must remain less than 70%). It is recommended that all optical accessories should be kept in a dry box when they are not being used. It is strongly recommended that a dehumidifier and air conditioning should be used continuously if the environment is too humid.

Cleaning: If there is any dust present on the optical surface, try to remove it with a fan or compressed air.

There are two recommended possible methods for removing fingerprints, grease or dust stains that cannot be removed by blowing with air:

- Blow gently onto the glass surface and then wipe with a clean cloth, paper for lenses or a cotton swab. Note that small cotton fibres may remain on the surface of the lens if a cotton swab is used.
- Use a cotton swab to clean the surface of the lens thoroughly. Do not use aggressive solvents.

The shelf life of the LED is of 5,000 hours. In case of non-function of the lighting and after checking that the batteries are charged, send the stereomicroscope back to Jeulin.

6. Cleaning and maintenance

The microscope must not be moved if it can be avoided. If it must be moved, the user must check that the eyepieces are well fixed into the eyepiece tubes before moving the instrument.

7. After-sales service

The guarantee is 2 years. The equipment must be returned to our workshops.

If you need any repairs, adjustments or spare parts, please contact:

JEULIN - Technical support
468 Rue Jacques-Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries